

TITLE

: URINE - THE ANALYSIS OF A SELECTION OF GROWTH PROMOTERS- LC-MS/MS



1 OBJECTIVE AND SCOPE

This SOP describes the analysis with LC-MS/MS of a selection of growth-promoting compounds: 17 α -trenbolone, 17 β -trenbolone, androsta-1,4-diene-3,17-dione (ADD), 17 β -boldenone, 17 α -boldenone, isoflupredone, dexamethasone, betamethasone, flumethasone, triamcinolone-acetonide, clobetasol, α -zeranol (α -zearalanol), β -zeranol (β -zearalanol), zearalanone, zearalenone, α -zearalenol, β -zearalenol, and 16 β -OH-stanozolol.

The method has a measuring range of 0 to 5.0 μ g/l.

2 DEFINITION

| | |
|----------------------------|---|
| MMS: | Matrix matched standards |
| ES | Electrospray |
| MS | Mass spectrometer |
| UPLC | Ultra Performance Liquid Chromatography |
| SPE | Solid Phase Extraction |
| ADD: | Androsta-1,4-diene-3,17-dione |
| 16 β -OH-stanozolol: | Metabolite of stanozolol |
| ACN: | Acetonitrile |
| MeOH | Methanol |
| HCOOH | Formic acid |
| Water | MilliQ H ₂ O |
| DMSO | Dimethyl sulfoxide |
| HAc | Acetic acid |
| NH ₃ | Ammonia |
| ISS | Internal standard solution |
| MSS | Mix standard solution |

3 PRINCIPLE

Glucuronated compounds in urine are deconjugated by glucuronidase at pH 6.5 – 7.5 and 37 °C. After clean-up with a SPE the sample is measured by means of LC-MS/MS.

4 CHEMICALS AND REAGENTS

Any reference to type and/or product is only to inform the user and identify the equipment and does not imply exclusion of similar equipment.

All reagents and chemicals must be at least pro analysis quality.

With 'water' is meant water, purified with a MilliQ® system with a minimum resistance of at least 18.2 M Ω .cm-1.

Instead of listed below, chemical products from other suppliers may be used, as long as they are from comparable quality.

4.1 Reference Standards

4.1.1 17 α -trenbolone (NMI D708)

4.1.2 17 β -trenbolone (Steraloids E3170-000)

4.1.3 17 β -Boldenone (Steraloids A200)

4.1.4 17 α -Boldenone (Rikilt, EU/CRL 53)

This is an ampoule containing 0.1 mg 17 α -boldenone

4.1.5 Isoflupredone (Research Plus 3005-16)

4.1.6 Flumethasone (Sigma F-9507)

4.1.7 Triamcinolone-acetonide (Sigma T-6501)

4.1.8 Betamethasone (Steraloids P0520-000)

4.1.9 Dexamethasone (Sigma D1756)

4.1.10 Clobetasol (Glaxo Wellcome)

4.1.11 ADD (androsta-1,4-diene-3,17-dione) (Steraloids A100)

4.1.12 α -Zeranol (Sigma Z-0292)

4.1.13 β -Zeranol (Sigma Z-0417)



European
Union
Reference
Laboratory

4.1.14 Zearalanone (TRC Z270450)

This is a vial containing 1 mg zearalanone.

4.1.15 Zearalenone (Sigma Z-2125)

4.1.16 α -zearalenol (Sigma Z0166)

4.1.17 β -zearalenol (Sigma Z 2000))

4.1.18 16β -OH-stanozolol NMI D621

This is an ampul containing 1 mg 16β -OH-stanozolol

4.1.19. Isoflupredone-d3 (d2 major) (TRC I816602)

This is a vial containing 1 mg Isoflupredone-d3 (d2 major).

4.1.20 17α -trenbolone-d3 (EURL 0101)

This is an ampoule containing 0.1 mg 17α -boldenone-d3.

4.1.21 17β -trenbolone-d3 (EURL 0102).

This is an ampoule containing 0.1 mg 17β -trenbolone-d3.

4.1.22 17β -Boldenone-d3 (EURL 0065)

This is an ampoule containing 0.1 mg 17β -boldenone-d3.

4.1.23 α -zearalanol-d4/ β -zearalanol-d4 (EURL 0006)

This is an ampoule containing 0.05 mg α -zearalanol-d4 and 0.05 mg / β -zearalanol-d4

4.1.24 α -zearalenol-d4 (EURL 0048)

This is an ampoule containing 0,1 mg α -zearalenol-d4.

4.1.25 β -zearalenol-d4 (EURL 0049)


This is an ampoule containing 0,1 mg β -zearalenol-d4

4.1.26 16β -OH-stanozolol-d3 (EURL 0064)

This is an ampoule containing 0.1 mg 16β -OH-stanozolol-d3.

4.1.27 Triamcinolone-acetonide-d6 (EURL 0063)

This is an ampoule containing 0.1 mg triamcinolone-acetonide-d6.



4.1.28 Dexamethasone-d4 (Cachesyn, CSTD361)

4.1.29 Zearalenone-d6 (TRC Z270502)

This is an vial containing 1 mg zearalenone-d6.

4.1.30 Zearalanone-d6 (TRC 270442)

This is an vial containing 0.25 mg zearalanone-d6.

4.1.31 17 β -Nortestosterone-D3

Nortestosterone-D3 (EURL 0004) 19-nortestosterone (16,16,17-D3)

This is an ampoule containing 0.1 mg nortestosterone-D3

4.2 Chemicals

4.2.1 Natriumhydroxide (Merck 1.06498)

4.2.2 Acetonitrile (Biosolve 01203502)

4.2.3 Acetic acid (Merck 1.00063)

4.2.4 β -Glucuronidase E Coli K12 (Roche 03707598001)

4.2.5 Methanol (Biosolv 13683502)

4.2.6 Disodium hydrogen phosphate dihydrate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$)(Merck 1.06580)

4.2.7 Potassium dihydrogen phosphate (KH_2PO_4)(Merck 1.04873)

4.2.8 Ammonium formate (Sigma 156264)

4.2.9 Formic acid (Merck 264)

4.2.10 Ethanol (Merck 1.00983)

4.2.11 Ammonia 25% (Merck 1.05432)

4.3 Standard solutions

Prepare standard solutions of individual compounds and internal standards. The concentrations, solvents, shelf lives (in freezer) and the amounts to weigh are given in table 1.



European
Union
Reference
Laboratory

Table 1: Solvent and concentrations of stock solutions

| Compound | Amount standard (mg) | Concentration stock solution (mg/l) | Solvent stock solution | Shelve lives (months) |
|--|----------------------|-------------------------------------|------------------------|-----------------------|
| 17 α -Trenbolone (4.1.1) | 2 - 10 | 1000 | methanol | 12 |
| 17 β -Trenbolone (4.1.2) | 2 - 10 | 1000 | methanol | 12 |
| 17 β -Boldenone (4.1.3) | 2 - 10 | 1000 | methanol | 12 |
| 17 α -Boldenone (4.1.4) | Ampul* | 10 | methanol | 12 |
| Isoflupredone (4.1.5) | 2 - 10 | 1000 | methanol | 12 |
| Flumethasone (4.1.6) | 2 - 10 | 1000 | methanol | 12 |
| Triamcinolone-acetonide (4.1.7) | 2 - 10 | 1000 | methanol | 12 |
| Betamethasone (4.1.8) | 2 - 10 | 1000 | methanol | 12 |
| Dexamethasone (4.1.9) | 2 - 10 | 1000 | methanol | 24 |
| Clobetasol (4.1.10) | 2 - 10 | 1000 | methanol | 12 |
| ADD (4.1.11) | 2 - 10 | 1000 | methanol | 12 |
| α -Zearalanol (4.1.12) | 2 - 10 | 1000 | methanol | 24 |
| β -Zearalanol (4.1.13) | 2 - 10 | 1000 | methanol | 24 |
| Zearalanone (4.1.14) | 1 | 100 | methanol | 24 |
| Zearalenone (4.1.15) | 2 - 10 | 1000 | methanol | 24 |
| α -Zearalenol (4.1.16) | 2 - 10 | 1000 | methanol | 24 |
| β -Zearalenol (4.1.17) | 2 - 10 | 1000 | methanol | 24 |
| 16 β -OH-stanozolol (4.1.18) | Ampul* | 100 | methanol | 24 |
| Isoflupredone-d3 (4.1.19) | 1 | 100 | methanol | 24 |
| 17 α -Trenbolone-d3 (4.1.20) | Ampul* | 10 | methanol | 24 |
| 17 β -Trenbolone-d3 (4.1.21) | Ampul* | 10 | methanol | 24 |
| 17 β -Boldenone-d3 (4.2.22) | Ampul* | 10 | ethanol | 24 |
| α -Zearalanol-d4/ β -zearalanol-d4 (4.1.23) | Ampul** | 10 | methanol | 24 |
| α -Zearalenol-d4 (4.1.24) | Ampul* | 10 | ethanol | 24 |
| β -Zearalenol-d4 (4.1.25) | Ampul* | 10 | ethanol | 24 |
| 16 β -OH-stanozolol-d3 (4.1.26) | Ampul* | 10 | methanol | 24 |
| Triamcinolone-acetonide-d6 (4.1.27) | Ampul* | 10 | methanol | 24 |
| Dexamethasone-d4 (4.1.28) | 2-10 | 1000 | methanol | 24 |
| Zearalenone-d6 (4.1.29) | 1 | 100 | methanol | 24 |
| Zearalanone-d6 (4.1.30) | 0.25 | 25 | methanol | 24 |
| 17 β -Nortestosterone-D3 (4.1.31) | Ampul* | 10 | methanol | 24 |

* Open the ampoule and add 1 ml solvent. Place during 1 minute in a ultra sonic bath, mix and transfer the contents by means of the solvent quantitative in a 10 ml flask. Fill up with solvent and mix.

** Open the ampoule and add 1 ml solvent. Place during 1 minute in a ultra sonic bath, mix and transfer the contents by means of the solvent quantitative in a 5 ml flask. Fill up with solvent and mix.

4.3.2 MSS I 10 mg/l

Pipet the amount of 100 μ l of each stock solution of 1000 mg/l and 1000 μ l of each stock of 100 mg/l in volumetric flasks of 10 ml and make up to volume with MeOH. The standard solutions are stored in a freezer for a maximum period of 12 months.

4.3.3 MSS II 0.1 mg/l

Pipet 100 μ l of MSS I (4.3.2), 100 μ l of each stock solution 10 mg/l and 10 μ l of each stock solution 100 mg/l in a volumetric flasks of 10 ml and make up to volume with MeOH. The standard solution is stored in a freezer for a maximum period of 3 months.

4.3.4 MSS III 0.02 mg/l

Pipet 2 ml of MSS II (4.3.3) in a volumetric flasks of 10 ml and make up to volume with MeOH. The standard solution is stored in a freezer for a maximum period of 3 months.

4.3.5 MSS IV 0.002 mg/l

Pipet 200 μ l of MSS II (4.3.3) in a volumetric flasks of 10 ml and make up to volume with MeOH. The standard solution is stored in a freezer for a maximum period of 3 months.

4.3.6 ISS dexamethasone-d4 10 mg/l

Pipet 100 μ l of stock solution 1000 mg/l in a volumetric flasks of 10 ml and make up to volume with MeOH. The standard solution is stored in a freezer for a maximum period of 12 months.

4.3.7 ISS I 0.1 mg/l

Pipet 100 μ l of ISS dexamethasone-d4 10 mg/l (4.3.6), 100 μ l of all other internal standards 10 mg/l, 40 μ l of internal standard 25 mg/l and 10 μ l of internal standard 100 mg/l in a volumetric flasks of 10 ml and make up to volume with MeOH. The standard solution is stored in a freezer for a maximum period of 12 months.

4.3.8 Working solution standards

Pipet 60 μ l of MSS II 0.1 mg/l (4.3.3) and 60 μ l of ISS I 0.1 mg/l (4.3.7) in a glass tube. Evaporate at 40 °C \pm 10% under a gentle stream of nitrogen until dry and re-dissolve in 120 μ l DMSO and 180 μ l 10% MeOH.

4.4 Reagents

4.4.1 HAc 8 N

Dilute 46 ml HAc to 100 ml with water. This solution is stable for 3 months at room temperature.

4.4.2 NaOH 8 N

Dissolve 32.0 gram NaOH in water and add water to a final volume of 100 ml. This solution is stable for 3 months at room temperature

4.4.3 Na₂HPO₄ 0.2 M

Dissolve 3.56 g Na₂HPO₄·2H₂O in water and add water to a final volume of 100 ml. This solution is stable for 1 month at room temperature.

4.4.4 KH₂PO₄ 0.2 M

Dissolve 2.72 g 4 KH₂PO₄ in water and add water to a final volume of 100 ml. This solution is stable for 1 month at room temperature.

4.4.5 0.2 M Phosphate buffer pH 7.0

Add KH₂PO₄ 0.2 M under stirring to Na₂HPO₄ 0.2 M till pH is 7.0 ± 0.2 (pH meter). This solution is stable for 1 month at room temperature.

4.4.6 Ammonium formate 1 M

Dissolve 6.3 g ammonium formate in water and add water to a final volume of 100 ml. This solution is stable for 3 months at room temperature

4.4.7 Wash solvent 1: 60% MeOH/2% HAc

Mix 38 ml water, 60 ml MeOH and 2 ml HAc. This solution is stable for 3 months at room temperature.

4.4.8 Wash solvent 2: 20% ACN/2% HAc

Mix 78 ml water, 20 ml ACN and 2 ml HAc. This solution is stable for 3 months at room temperature.

4.4.9 Wash solvent 3: 20% ACN

Mix 80 ml water with 20 ml ACN. This solution is stable for 3 months at room temperature.

4.4.10 Wash solvent 4: 10% ACN/2% NH₃

Mix 82 ml water, 8 ml ammonia (25%) and 10 ml ACN. This solution is stable for 3 months at room temperature.

4.4.11 Wash solvent 5: 50% MeOH/2% NH₃

Mix 42 ml water, 8 ml ammonia (25%) and 50 ml MeOH. This solution is stable for 3 months at room temperature.

4.4.12 Wash solvent 6: 60% MeOH

Mix 40 ml water and 60 ml MeOH. This solution is stable for 3 months at room temperature.

4.4.13 Mobile phase A: water/ACN/HCOOH/ammonium formate 1 M 900/100/0.02/2 (v/v/v/v)

Mix 900 ml water, 100 ml ACN, 0.02 ml HCOOH and 2 ml ammonium formate 1 M (4.4.6). This solution is stable for 1 month at room temperature.

4.4.14 Mobile phase B: water/ACN/HCOOH/ammonium formate 1 M 100/900/0.02/2 (v/v/v/v)

Mix 100 ml water, 900 ml ACN, 0.02 ml HCOOH and 2 ml ammonium formate 1 M (4.4.6). This solution is stable for 1 month at room temperature.

5 EQUIPMENT

Any reference to type and/or product is only to inform the user and identify the equipment and does not imply exclusion of similar equipment.

5.1 Analytical balance with an accuracy of 0.02 mg or better (Mettler AT 261)

5.2 pH indicator strips 6.5 – 10,0 (Merck 9543)

5.3 Water bath adjustable at 37 °C (Julabo MD/25B)

5.4 Glass tube 12 ml (Beldico, 8739007)

5.5 Polypropylene tube 12 ml with screwcap (Greiner 163275)

5.6 Evaporation unit (Turbo-Vap LV Zymark 44467)

5.7 Evaporator for 96 wells (Caliper LifeSciences Turbovap 96)

5.8 Vortex (IKA Vibrofix VF1)

5.9 Vortex (IKA Genius 3)

5.10 OASIS HLB LP 96-well plate 60 µm (60mg) (Waters 186000679)

5.11 Extraction plate manifold for Oasis 96-wells plates (Waters 186001831)

5.12 Vacuum pump (KNF Lab Laboport D79112)

5.13 96-Well collection plates for Waters Acquity UPLC I-Class, well volume 2 ml (WAT 058985)

5.14 96-Well PTFE/Silicone seal with pre-slit (Waters 186006335)

5.15 12-Channel pipette (Rainin 100-1200 µl LTS Pipet Lite, L12-1200XLS, Mettler Toledo 17014497)

5.16 pH meter (Scott CG 840)

5.17 LC-MS opstelling

5.17.1 Waters Acquity pump and injection system

5.17.2 Analytical kolom: Acquity UPLC BEH C₁₈, 100 mm x 1.0 mm, 1.7 µm (Waters 186002346)

5.17.3 Mass spectrometer: Waters, Xevo TQS with ESI interface.

6 PROCEDURE

6.1 General

This SOP describes the quantification and confirmation of a selection of growth promoters in cattle urine.

6.2 Safety precautions

Take precautions to prevent inhalation and contact with skin from standards and chemicals. Work in a fume hood and wear gloves if necessary (SOP F0071). [2].

6.3 Sample pre-treatment

Frozen samples are thawed at room temperature or overnight in the fridge. Homogenize samples before a test unit is taken. After that place sample as soon as possible back in the freezer.

6.4 Test unit

Take 0.5 ml of urine.

6.5 Description procedure

A checklist of the method is given in appendix 2.

6.5.1 Composition of samples series

A serie consists of samples, 7 samples with spike before processing (MMS) and a 1st-line control. For the MMS, and 1st-line control a urine sample is used which in previous analysis no compounds were detected. For the 1st-line control another blank urine is taken than used for MMS.

6.5.2 Preparing MMS, 1st-line control and urine samples

6.5.2.1 Preparing MMS

Pipet in separate polypropylene tubes the MSS and ISS according to table 2.

Table 2. Preparing MMS

| Matrix Matched Standards | μl MSS IV (0.002 mg/l)(4.3.5) | μl MSS III (0.02 mg/l)(4.3.4) | μl MSS II (0.1 mg/l)(4.3.3) | μl ISS I 0.1 mg/l (4.3.7) |
|------------------------------|---|---|---|---|
| MMS A (0 $\mu\text{g/l}$) | 0 | | | 25 |
| MMS B (0.1 $\mu\text{g/l}$) | 25 | | | 25 |
| MMS C (0.2 $\mu\text{g/l}$) | 50 | | | 25 |
| MMS D (0.5 $\mu\text{g/l}$) | | 12.5 | | 25 |
| MMS E (1.0 $\mu\text{g/l}$) | | 25 | | 25 |
| MMS F (2.0 $\mu\text{g/l}$) | | 50 | | 25 |
| MMS G (5.0 $\mu\text{g/l}$) | | | 25 | 25 |

Evaporate under a gently stream of nitrogen at $40\text{ }^{\circ}\text{C} \pm 10\%$ until just dry. Pipet in every tube a test unit of a blank sample. Mix with vortex. Proceed from 6.5.3.

6.5.2.2 1st-line control

Pipet in polypropylene tube 25 μl of ISS I (0.1 mg/l)(4.3.7) and 25 μl MSS III (0.02 mg/l)(4.3.4). Evaporate under a gently stream of nitrogen at $40\text{ }^{\circ}\text{C} \pm 10\%$ until just dry. Pipet 1 test unit of another blank urine than used for MMS in the polypropylene tube. Mix with vortex. Proceed from 6.5.3.

6.5.2.3 Preparing samples

Pipet for each sample in separately polypropylene tubes 25 μl of ISS I (0.1 mg/l)(4.3.7). Evaporate under a gently stream of nitrogen at $40\text{ }^{\circ}\text{C} \pm 10\%$ until just dry. Pipet for each samples a test unit in the polypropylene tube. Mix with vortex. Proceed from 6.5.3.

6.5.3 Deconjugation

Add to each tube 500 μl 0.2 M phosphate buffer pH 7.0 (4.4.5) and vortex. Check the pH with pH indicator strips and adjust if necessary with 8M NaOH and/or 8M HAC or diluted solutions thereof until the pH 6.5-7.5 is. Add 5 μl β -glucuronidase and mix. Place in a water bath at $37\text{ }^{\circ}\text{C}$ during 16 hours. Cool to room temperature.

6.5.4 Sample clean-up by 96-wells SPE

Conditioning a SPE plate with 1 ml of MeOH followed by 1 ml of water. Bring whole sample in SPE well. Wash with 1 ml water and dry. Wash successively with 1 ml wash solvent 1 to 6. Dry after each wash step. Elute with 1 ml ACN. Collect eluate in a 96-well collection plate which in each well is already 20 μl of DMSO. Place plate in a 96-wells evaporator and dry at $60\text{ }^{\circ}\text{C} \pm 10\%$ till only DMSO is left. Add 30 μl of 10% MeOH and vortex (5.9) carefully. The samples are now suitable for LC-MS/MS analysis.

6.5.5 LC-MSMS analysis

6.5.5.1 LC-conditions

| | |
|-----------------------|---|
| Column : | UPLC BEH C ₁₈ , 100 mm x 1.0 mm, 1.7 μm |
| Mobile phase A: | water/ACN/HCOOH/ammonium formate 1 M 900/100/0.02/2 (v/v/v/v) |
| Mobile phase B: | water/ACN/HCOOH/ammonium formate 1 M 100/900/0.02/2 (v/v/v/v) |
| Column temperature: | 60 °C |
| Vialtray temperature: | 20°C |
| Injection volume | 5 μl |
| Flow: | 0.15 ml/min |
| Sample manager purge: | H ₂ O/ACN 700/300 (v/v) |
| Sample manager wash: | ACN/water/HCOOH 900/100/0.1 (v/v/v) |
| Seal wash: | H ₂ O/MeOH 900/100 (v/v) |
| Gradient: | see table 3 |

Table 3: Gradient

| Time (min) | Mobiele fase A (%) | Mobiele fase B (%) | Flow (ml/min) |
|------------|--------------------|--------------------|---------------|
| 0.0 | 80 | 20 | 0.15 |
| 0.2 | 80 | 20 | 0.15 |
| 3.2 | 70 | 30 | 0.15 |
| 7.5 | 50 | 50 | 0.15 |
| 7.6 | 0 | 100 | 0.15 |
| 8.6 | 0 | 100 | 0.15 |
| 8.7 | 80 | 20 | 0.15 |
| 10.0 | 80 | 20 | 0.15 |

6.5.5.2 MS-conditions

The MS parameters are guidelines. To obtain a better performance these parameters can be changed.

Ionisation mode: ESI, negative and positive

| Parameter | ESI, positive | ESI, negative* |
|-------------------|---------------|----------------|
| Capillair voltage | 3.0 | 2.2 |
| LM 1 resolution | 3.0 | 2.8 |
| HM 1 resolution | 15.0 | 15.0 |
| Ion energy 1 | 1.0 | 1.9 |
| LM 2 resolution | 3.0 | 3.0 |
| HM 2 resolution | 15.0 | 15.0 |
| Ion energy 2 | 1.0 | 1.3 |

* = specific values for TQS1

| | |
|--------------------------|---|
| Cone: | compound depends V |
| Source temperature: | 150 °C |
| Desolvation temperature: | 400 °C |
| Desolvation gas flow: | 800 l/hour |
| Cone gas flow: | 150 l/hour |
| Nebuliser gas flow | 7 bar |
| CID gas: | argon, p = 2,2 * 10 ⁻³ mbar, purity > 99.998% |
| CID gas flow: | 0.18 ml/min |

The compounds fragment to structure related product ions. In table 4 the theoretical mono-isotopic masses of the precursor ions and corresponding product ions are displayed. For each of the mono-isotopic masses the permitted deviation is ± 0.5 Da.

Table 4: Monitored MRM's

| Component | Precursor ion (m/z) | Product ion (m/z) | Cone Voltage (V) | Collision energy (eV) | ESI |
|-------------------------------------|---------------------|-------------------|------------------|-----------------------|-----|
| Isoflupredone* | 423.2 | 293.0 | 20 | 30 | - |
| | | 347.1 | | 15 | |
| Isoflupredone-d3* | 426.2 | 350.1 | 20 | 17 | - |
| Dexamethasone/Betamethasone* | 437.2 | 307.0 | 20 | 30 | - |
| | | 361.2 | | 15 | |
| Dexamethasone-d4* | 441.2 | 363.10 | 20 | 15 | - |
| Flumethasone* | 455.2 | 305.0 | 20 | 30 | - |
| | | 379.1 | | 20 | |
| 17β-Trenbolone/17α-Trenbolone | 271.2 | 199.1 | 30 | 25 | + |
| | | 253.2 | | 20 | |
| 17β-Trenbolone-d3/17α-Trenbolone-d3 | 274.1 | 256.2 | 30 | 20 | + |
| 17β-Boldenone/17α-Boldenone | 287.2 | 121.1 | 25 | 25 | + |
| | | 135.1 | 30 | 15 | |
| 17β-Boldenone-d3 | 290.2 | 121.1 | 25 | 25 | + |
| ADD | 285.2 | 121.1 | 25 | 25 | + |
| | | 147.1 | 15 | 15 | |

| Component | Precursor ion (m/z) | Product ion (m/z) | Cone Voltage (V) | Collision energy (eV) | ESI |
|--|---------------------|-------------------|------------------|-----------------------|-----|
| 16 β -OH-Stanozolol | 345.2 | 81.1 | 30 | 45 | + |
| | | 95.1 | | 40 | |
| 16 β -OH-Stanozolol-d3 | 348.3 | 81.1 | 50 | 45 | + |
| Triamcinolone-acetonide | 435.3 | 339.2 | 40 | 13 | + |
| | | 415.3 | | 8 | |
| Triamcinolone-acetonide-d6 | 441.3 | 421.3 | 40 | 8 | + |
| Clobetasol | 411.2 | 355.1 | 20 | 13 | + |
| | | 391.2 | | 5 | |
| Zearalenone | 317.2 | 131.2 | 44 | 30 | - |
| | | 175.1 | | 24 | |
| β -Zearalenol/ α -Zearalenol/Zearalanone | 319.2 | 160.2 | 30 | 30 | - |
| | | 275.2 | | 21 | |
| β -Zearalenol-d4/ α -zearalenol-d4 | 323.2 | 160.1 | 45 | 32 | - |
| Zearalanone | 319.2 | 205.2 | 54 | 22 | |
| | | | | | |
| β -Zearalanol/ α -Zearalanol | 321.2 | 259.2 | 56 | 26 | - |
| | | 277.2 | | 20 | |
| β -Zearalanol-d4/ α -Zearalanol-d4/Zearalanone-d6 | 325.1 | 281.1 | 50 | 20 | - |
| 17 β -Nortestosterone-d3 | 278.2 | 109.1 | 30 | 25 | + |

* = formic acid adduct

6.5.5.3 Initial test LC system

By means of working standard solution (4.3.8) the retention time and stability of the system are tested. Inject hereafter MMS B. Calculate the signal/noise ratio of the lowest ion of each compound. The signal/noise ratio should be at least 6 for each compound.

6.5.5.4 Injection order

The following injection order is advisable:

- MMS B (initial test)
- blank
- MMS A - G
- blank
- 1st-line control
- blank
- samples

-blank

- MMS A – G

7 RESULTS

7.1 Calculations

Table 5 is an overview of which internal standard is used for calculations.

Table 5: Compound with internal standard used for calculations

| Compound | Internal standard |
|-----------------------------|----------------------------|
| ADD | 17β-Nortestosterone-D3 |
| Isoflupredone | Isoflupredone-d3 |
| Dexamethasone/Betamethasone | Dexamethasone-d4 |
| Flumethasone | Dexamethasone-d4 |
| 17β-Trenbolone | 17β-Trenbolone-d3 |
| 17α-Trenbolone | 17α-Trenbolone-d3 |
| 17β-Boldenone | 17β-Boldenone-d3 |
| 17α-Boldenone | 17β-Boldenone-d3 |
| Triamcinolone-acetonide | Triamcinolone-acetonide-d6 |
| 16β-OH-Stanozolol | 16β-OH-stanozolol-d3 |
| Clobetasol | Dexamethasone-d4 |
| Zearalenone | α-Zearalenone-d6 |
| β-Zearalenol | β-Zearalenol-d4 |
| α-Zearalenol | α-Zearalenol-d4 |
| Zearalanone | α-Zearalanone-d6 |
| β-Zearalanol | β-Zearalanol-d4 |
| α-Zearalanol | α-Zearalanol-d4 |

Equation I: Calculation response factor (RF)

$$RF = \left(\frac{Area_{compound}}{Area_{IS}} \right)$$

With:

RF = response factor

$Area_{compound}$ = sum area product ions

$Area_{IS}$ = area product ion internal standard

Equation II: Calculation of the concentration in the sample (X)

$$X = \left(\frac{RF - b}{a} \right)$$

With:

X = concentration of a compound in the sample ($\mu\text{g/l}$)

RF = response factor (equation I (7.1)).

b = intersection point of the calibration curve with the y-axis (follows from linear regression*)

a = slope of the calibration line (follows from linear regression*)

*Plot the area of the MMS measured before and after the sample extracts as function of the added levels. Calculate the linear regression using the least squares method.

Equation III: Calculation relative retention time (RRT)

$$RRT = \left(\frac{RT_{\text{compound}}}{RT_{IS}} \right)$$

With:

RRT = relative retention time

RT_{compound} = retention time of compound

RT_{IS} = retention time of internal standard

Equation IV: Calculation of the relative deviation of the relative retention time (ΔRRT)

$$\Delta RRT = \left(\frac{RRT_{\text{sample}} - RRT_{\text{mean}}}{RRT_{\text{mean}}} \right) \times 100\%$$

With:

ΔRRT = relative deviation of the relative retention time of compound compare to mean retention time of MMS B up to and including G before and after samples (%)

RRT_{sample} = relative retention time of a compound in the sample

RRT_{mean} = mean relative retention time of a compound in the MMS B up to and including G before and after samples

Equation V: Calculation minimal respons of internal standard for assurance determination

$$A_{\text{min}} = \frac{6 * A_{\text{MMS_RL}}}{S / N}$$

With:

A_{min} = minimal of internal standard in a sample where level of $\frac{1}{4}$ of the proposed target is guaranteed

$A_{\text{MMS_RL}}$ = lowest area of internal standard of MMS B before or after samples

S/N = signal/noise ratio of the lowest product ion of MMS B before the samples

Equation VI: Ion ratio (R):

$$R = \left(\frac{A_{low}}{A_{high}} \right) \times 100\%$$

With:

- R = ion ratio (%)
- A_{low} = area of the product ion with the lowest intensity
- A_{high} = area of the product ion with the highest intensity

Equation VII: Relative deviation of the ion ratio (D):

$$\Delta R = \left(\frac{R_{sample} - R_{mean.}}{R_{mean}} \right) \times 100\%$$

With:

- D = relative deviation of the ion ratio of the compound in the sample compared to the average ion ratio of the compound in MMS B up to and including G before and after samples
- R_{sample} = ion ratio of the compound in the sample
- R_{mean.} = mean ion ratio of the compound in MMS B up to and including G before and after samples

7.2 Criteria

7.2.1 General

For acceptance of the analysis results, the performance criteria (stated below) have to be met. The concentrations are automatically calculated when the result sheet (appendix 1) is used. At the same time the excel sheet will show if the calculated analysis results meet the stated criteria.

7.2.2 Drift in sensitivity

The MMS analysed before the sample extracts compared with the MMS analysed after the sample extracts are indicative for the performance of the LC-MS/MS system. For both MMS series, the sum of the area of the product ions is plotted as function of the added concentration to the sample (µg/l). Two slopes for both MMS series are constructed, which may not differ by more than 25%.

7.2.3 Linearity

The MMS are used to determine the linearity of the LC-MS/MS system and to determine if the sample pre-treatment is done correctly. For both MMS series, the sum of the area of the product ions is plotted as function of the added concentration in the sample (µg/l). Apply linear regression using the least squares method. The correlation of both lines should be ≥ 0.990.

7.2.4 Sensitivity

The S/N of the least intense product ion of MMS 0.25*RW should be ≥ 6 .

7.2.5 Accuracy

Calculate the concentration in the 1st line control using equation II. Calculate the accuracy for each compound by dividing the concentration by the actual spiked level (%). Fill in the Shewart charts. SOP F0057 discusses Shewhart charts.

7.2.6 Maximum deviation retention time

Determine the average retention time of the compounds in MMS B - G before and after samples the sample extracts. The deviation in the individual retention times may not differ more than 2.5% compared to the mean retention time of the compounds in the MMS B - G.

7.2.7 Maximum deviation ion ratio

Calculate the mean ion ratio of the MMS B up to and including G before and after the samples. The deviation in the ion ratio of each MMS sample compared to the mean should not exceed the EU-criteria in table 6.

Table 6: Maximum permitted tolerances for relative ion intensities according to EU criteria [3]

| Average ion ratio MMS (R) | Permitted relative deviation (D) |
|------------------------------|-------------------------------------|
| R > 50% | $\leq 20\%$ |
| 20% < R \leq 50 % | $\leq 25\%$ |
| 10% < R \leq 20 % | $\leq 30\%$ |
| R \leq 10% | $\leq 50\%$ |

7.2.8 Quantification

The concentration of the compounds in the sample is expressed as $\mu\text{g/l}$ and is calculated using the linear regression curves calculated from the MMS series analysed before and after the sample extracts. Calculate the concentrations using equation I and II (7.1). The concentration is automatically calculated using the result sheet (appendix 1).

7.2.9 Criteria for individual samples

Calculate for each sample the relative deviation of the retention time and the relative deviation of the ion ratio. When the relative deviation of the retention time does not exceed 2.5% and the relative deviation of the ion ratio does not exceed the EU-criteria and the concentration exceeds $CC\alpha$, the sample is non-compliant. The second and/or third validation level decides whether the sample is reanalyzed.

7.2.10 Accuracy

Calculate the concentration of a compound in the 1st-line control sample using equation II. Calculate the accuracy for each compound by dividing the concentration by the actual spiked level (%). Fill in the Shewart charts. SOP F0057 discusses Shewart charts

7.3 Final results


When a compound is identified in the sample with a concentration exceeding $CC\alpha$, the concentration is reported in LIMS. When no compound is identified in the sample "<0.25 RW" (lowest MMS) is reported in LIMS.

8 REGISTRATION

All information concerning the analysis is described in the labjournal. The MS/MS data are digital archived.

LITERATURE

- [1] EU Beslissing 2002/657/EC, J.EUR. Communities, L221 (2002) 8
- [2] SOP F0057 – Controle kaarten – Gebruik en beheer
- [3] SOP F0071 – Persoonlijke beschermingsmiddelen binnen het RIKILT



Appendix 1 belonging to SOP A1185

BULAGE bij SOP A1185

Waarnemingsformulier voor bepaling en confirmatie versie 1

Analysedatum: _____
 Matrie: **Mine**
 Component: **Zwavelsteen** rapportagepers. (dd-mm-jjjj): _____

Initiële test

| Data file | Omschrijving | RT | Area m/z hoogste intensiteit | Area m/z laagste intensiteit | Concentratie laag | Concentratie hoog |
|-----------|--------------|----|------------------------------|------------------------------|-------------------|-------------------|
| | | | | | 1.0 | 10.0 |

Criteria initiële test

| Omschrijving | S/N m/z laagste intensiteit | Criterium | Volgdat | Serie geaccepteerd | Paraat |
|----------------|-----------------------------|-----------|---------|--------------------|--------|
| Gemiddelde m/z | | S/N > 6 | | | |

MMS

MMS geanalyseerd vóór monsters

| Data file | Omschrijving | RT | Area m/z hoogste intensiteit | Area m/z laagste intensiteit | IS RT | Area IS | Concentratie laag (µg/L) | Concentratie hoog (µg/L) | Relatieve RT | Respons factor | Afwijking m/z (RT (%)) | Isotratio | Afwijking isotratio (%) | Gehalte (µg/kg (aq.)) | Waarde |
|-----------|--------------|----|------------------------------|------------------------------|-------|---------|--------------------------|--------------------------|--------------|----------------|------------------------|-----------|-------------------------|-----------------------|---------------|
| | | | | | | | 0.0 | 0.1 | | | | | | | Richting |
| | | | | | | | 0.2 | 0.2 | | | | | | | Stippunt y-as |
| | | | | | | | 0.5 | 0.5 | | | | | | | Correlatie |
| | | | | | | | 1.0 | 1.0 | | | | | | | |
| | | | | | | | 2.0 | 2.0 | | | | | | | |
| | | | | | | | 5.0 | 5.0 | | | | | | | |
| | | | | | | | 10.0 | 10.0 | | | | | | | |

Gemiddelde area IS vóór: _____

MMS geanalyseerd ná monsters

| Data file | Omschrijving | RT | Area m/z hoogste intensiteit | Area m/z laagste intensiteit | IS RT | Area IS | Concentratie laag (µg/L) | Concentratie hoog (µg/L) | Relatieve RT | Respons factor | Afwijking m/z (RT (%)) | Isotratio | Afwijking isotratio (%) | Gehalte (µg/kg (aq.)) | Waarde |
|-----------|--------------|----|------------------------------|------------------------------|-------|---------|--------------------------|--------------------------|--------------|----------------|------------------------|-----------|-------------------------|-----------------------|---------------|
| | | | | | | | 0.1 | 0.1 | | | | | | | Richting |
| | | | | | | | 0.2 | 0.2 | | | | | | | Stippunt y-as |
| | | | | | | | 0.5 | 0.5 | | | | | | | Correlatie |
| | | | | | | | 1.0 | 1.0 | | | | | | | |
| | | | | | | | 2.0 | 2.0 | | | | | | | |
| | | | | | | | 5.0 | 5.0 | | | | | | | |
| | | | | | | | 10.0 | 10.0 | | | | | | | |

Gemiddelde area IS ná: _____

Verloop gevoeligheid serie

| Verloop responsfactor | Gemiddelde area IS | Gem. rel. RT | Staprel. RT | RSD rel. RT | Gem. isotratio | Staprel. isotratio | RSD isotratio |
|-----------------------|--------------------|--------------|-------------|-------------|----------------|--------------------|---------------|
| | | | | | | | |

MMRS

| Data file | Omschrijving | RT | Area m/z hoogste intensiteit | Area m/z laagste intensiteit | IS RT | Area IS | Concentratie laag (µg/L) | Concentratie hoog (µg/L) | Relatieve RT | Respons factor | Afwijking m/z (RT (%)) | Isotratio | Afwijking isotratio (%) | Telvoetgeving (µg) |
|-----------|--------------|----|------------------------------|------------------------------|-------|---------|--------------------------|--------------------------|--------------|----------------|------------------------|-----------|-------------------------|--------------------|
| | | | | | | | 0.50 | | | | | | | |

Criteria voor runacceptatie:

| Omschrijving | Waarde | Criterium | Volgdat | Serie geaccepteerd | Paraat |
|------------------------------------|--------|-----------|---------|--------------------|---------------------|
| Verloop responsfactor | | ≤ 20.0% | | | 1e validatierversie |
| Lineairiteit | | ≥ 0.990 | | | Paraat |
| Gevoeligheid m/z laagste gevoeging | | SNR ≥ 6 | | | |
| Max. afw. isotratio MMS | | ≤ 2.5% | | | 2e validatierversie |
| Max. afw. m/z RT MMS | | ≤ 2.5% | | | Paraat |

Criteria voor acceptatie monster resultaat:

| Omschrijving | Criterium |
|---|-----------|
| Min. signaal IS voor bepaling detectie op rapportagegrens | |
| Max. afw. RT conform EU-criteria | ≤ 2.5% |
| Max. afw. isotratio conform EU-criteria | ≤ 2.5% |

Monsters

| Data file | Omschrijving | RT | Area m/z hoogste intensiteit | Area m/z laagste intensiteit | IS RT | Area IS | Bepaling det. 1-4 MMS | S/N Area m/z laagste intensiteit | Relatieve RT | Afwijking m/z (RT (%)) | Respons factor | Isotratio | Afwijking isotratio (%) | Gehalte (µg/kg (aq.)) | Confirmatie resultaat |
|-----------|--------------|----|------------------------------|------------------------------|-------|---------|-----------------------|----------------------------------|--------------|------------------------|----------------|-----------|-------------------------|-----------------------|-----------------------|
| | | | | | | | | | | | | | | | Bepaling det. |

Opmerkingen

Lineairiteitsplot

Date: _____
 Name analyst: _____
 Labjournal/page: _____
 Working list number(s) _____

- Preparing MMS. 1st-line control and samples

| Matrix Matched Standards | μl MSS IV (0.002 mg/l)(4.3.5) | μl MSS III (0.02 mg/l)(4.3.4) | μl MSS II 0.1 mg/l)(4.3.3) | μl ISS 0.1 mg/l (4.3.7) |
|------------------------------|---|---|--|---------------------------------------|
| MMS A (0 $\mu\text{g/l}$) | 0 | | | 25 |
| MMS B (0.1 $\mu\text{g/l}$) | 25 | | | 25 |
| MMS C (0.2 $\mu\text{g/l}$) | 50 | | | 25 |
| MMS D (0.5 $\mu\text{g/l}$) | | 12.5 | | 25 |
| MMS E (1.0 $\mu\text{g/l}$) | | 25 | | 25 |
| MMS F (2.0 $\mu\text{g/l}$) | | 50 | | 25 |
| MMS G (5.0 $\mu\text{g/l}$) | | | 25 | 25 |
| 1st-line control | | 25 | | 25 |
| Samples | | | | 25 |

- Evaporate until just dry at 40 °C
- Pipet 0.5 ml blank urine (1st line control another urine than MMS)
- Pipet 0.5 ml of each sample in tube
- Add 500 μl 0.2 M phosphatebuffer pH 7.0
- Mix and adjust to pH 6.5 - 7.5 with HAC/NaOH
- Add 5 μl β -Glucuronidase and mix
- Place during 16 hours in water bath at 37 °C
- Cool to room temperature
- Conditioning a SPE plate with 1 ml of MeOH followed by 1 ml of water
- Transfer whole sample in SPE well
- Wash with 1 ml of water and dry
- Wash with 1 ml wash solvent 1 and dry
- Wash with 1 ml wash solvent 2 and dry
- Wash with 1 ml wash solvent 3 and dry
- Wash with 1 ml wash solvent 4 and dry
- Wash with 1 ml wash solvent 5 and dry
- Wash with 1 ml wash solvent 6 and dry
- Pipet 20 μl DMSO in each well of the collection plate
- Elute with 1 ml ACN.
- Place plate in a 96-wells evaporator and dry at 60 °C \pm 10% till only DMSO is left.

- Add 30 μ l of 10% MeOH.and vortex (5.9) carefully
- The samples are now suitable for LC-MS/MS analysis

