

1. General information	1
1.1 pH stuff	1
1.2 Abbreviations and solutions	2
1.3 Tips and ways to reduce the amount of Keratins in your samples.....	2
1.4 Recommended procedures	3
1.5 Sample losses, methods comparison.....	3
1.6 nLC-MSMS sample necessities.	4
1.7 General sample cleanup procedures with μ Columns.	4
1.8 About Methionine oxidation	4
2. Protein determination (BCA).....	5
3. Gel free protein digestion methods.....	6
3.1 In-StageTip (iST) sample preparation for maximally 5 ug protein without detergent.	6
3.2 Filter aided sample preparation (FASP) for maximally 10 ug protein	6
3.3 Normal "In solution" trypsin digestion	7
3.4 Methanol and TriFluoroEthanol (TFE) sample preparation method.....	7
3.5 Other methods.....	7
4. In-Gel Digestion method (IGD)	8
4.1 General info	8
4.2 Recommended procedure for CCB or Oriole protein gel staining:.....	8
4.2.1 Colloidal Coomassie Preparing Staining Solution	8
4.2.2 Colloidal Coomassie Gel staining procedure	8
4.3 In Gel Digestion protocol	9
4.3.1 Remarks	9
4.3.2 Procedure	9
5. Quantitation	10
5.1 Relative quantitation (= Sample compared to Control)	10
5.1.1 Relative quantitation by on column peptide dimethyl labelling protocol.....	11
5.2 Absolute quantitation	12
6. Phosphopeptide enrichment methods (S and T only)	13
6.1 Phosphopeptides (S, T) enrichment by Titanium Dioxide (TiO ₂) Chromatography.....	14
6.1.1. Enrichment with normal TiO ₂ in μ Columns.....	14
6.1.2. Enrichment with magnetic Ti ⁴⁺ beads.....	14
6.2 Phosphopeptide (S,T) sample preparation by sequential elution from IMAC.	15



1. General information

Do not use any chemicals that react with amino acids for proteomics sample prep (except for the ones described in this protocol to alkylate cysteines). Avoid e.g. popular (irreversibly reacting) protease inhibitors like Aminoethylbenzenesulfon (AEBS, Y + 183.0354) or E-64 in some protease inhibitor cocktails and many phosphatase inhibitors. Instead use reversible inhibitors when necessary e.g.

http://www.sigmaaldrich.com/catalog/product/sigma/mssafe?lang=en®ion=NL&cm_sp=Insite--recent_fixed--recent5-2

Inhibitors info: <https://www.thermofisher.com/nl/en/home/life-science/protein-biology/protein-biology-learning-center/protein-biology-resource-library/pierce-protein-methods/protease-phosphatase-inhibitors.html>.

1.1 pH stuff

You can check the pH of your sample by putting 0.1 ul of sample (or less) on a piece of pH paper.

pH UP from pH 2 to pH 8

1 ml/l HCOOH = formic acid (FA): pH = 2.4

For 1 ml: Make pH 8 – 8.5 by adding 15 – 20 ul 10* diluted conc. NH₃ (max 37%).

0.5 ml/l TFA = TriFluoroAcetic acid: pH = 2.1

For 1 ml: Make pH 8 – 8.5 by adding 6- 9 ul 10* diluted conc. NH₃ (max 37%).

pH DOWN from pH 8 to pH 2 - 4

50 mM ABC pH 8.0:

For 1 ml: Make pH 3 by adding 35 ul 10* diluted conc. TFA.

1.2 Abbreviations and solutions

AcNi	acetonitril		
ABC	Ammonium BiCarbonate	NH_4HCO_3 (lab 2)	(M=79) 50 mM ABC pH 8.0 = 0.2 g / 50 ml
AmAc	Ammonium Acetate	$\text{NH}_4\text{CH}_3\text{COOH}$	(M = 77) 10 mM = 38.5 mg / 50 ml
TFE	TriFluoroEthanol (lab 2, safety cupboard MS stuff)		
TFA	TriFluoro-Acetic acid	! FUMING = below fume hood = add in fume cupboard only	
TCEP	reductor Tris(CarboxyEthyl)Phosphine(-20°C)	M=287	100 mM = 28.7 mg / ml
Urea	(always prepare fresh! Do not warm up!) (weighing room)	M=60	8 M = 0.48 g / ml
DTT = Dithiotreitol	(4°C) FRESH! in 50 mM ABC	M=154	50 mM = 7.7 mg/ml, 500 mM = 77 mg/ml
ME = β -MercaptoEthanol	(4°C) FRESH! Concentrated = 14.3 M	M=78.1	50 mM = 3.5 ul/ml
AA = AcrylAmide	(Wear gloves) FRESH! in 50 mM ABC	M=71.1	50 mM = 3.6 mg/ml
IAA = Iodoacetamide	(4°C) FRESH! in 100 mM Tris pH 8	M=185	100 mM = 19 mg/ml 500mM = 92.5 mg/ml
Cysteine (Fluka 30090, >99%, lab 2)	in 50 mM ABC	M=121	125 mM = 15 mg/ml 500 mM = 60 mg/ml

Stock Tris (10*):	1 M Tris bring to pH 8 with HCl	M=121	1 M = 121 g/l = 6.0 g / 50 ml
ST buffer stock:	4%(w/v) SDS in 100 mM Tris/HCl pH 7.6		100 ul 1M Tris + 0.04 g SDS made up to 1.0 ml
SDT-lysis buffer:	50 mM DTT in ST buffer	M=154	7.7 mg / 1 ml ST buffer
UT: 8 M urea (Sigma, U0631) in 100 mM Tris/HCl pH 8		M=60	100 ul 1M Tris + 0.48 g urea made up to 1.0 ml pH will increase to 8.2 due to the addition of urea.

Trypsin: we have a stock solution of Bovine Sequencing grade Trypsin (Roche 11 047 841 001) of 0.5 ug/ul (500 ng/ul) in 1 mM HCl. Generally dilute 100* in ABC before use.

- 0.5 ml protein LoBind tube: order# eppe0030108.094
- 2.0 ml protein LoBind tube: order# eppe0030108.132
- Pall 3K or 10K omega filter (20kDa / 50kDa cutoff) order# Pall OD003C34 = Sigma-Aldrich [Z722049-100EA](#)

1.3 Tips and ways to reduce the amount of Keratins in your samples.

- DO:** Use commercial electrophoresis gels. They contain less keratin than your own ones.
Use Eppendorf low binding tubes to minimize protein losses (see figure below).
Keep gels covered as much as possible. Put them into a new square petri dish for staining.
- Use Nitril gloves, not Latex ones.
- Wash your hands under running tap water before you start and as often as possible in between handlings.
- Do not work in a standard flow cabinet intended for microbiological work. They increase the amount of keratin passing your vials resulting in more keratin in your samples.
- Do not wear clothes of wool.
- Try not to lean over the samples too much.
- Never ever use glass vials for proteins. You will lose a lot protein.
- Use new (or only used for proteomics) throw away (polypropylene) plastics.
- Do not use hand creams when preparing samples for MS.
- Do not use any glassware that has been cleaned with detergent (e.g. in a washing machine).
- Cheap (non-Eppendorf) micro tubes may contain polymers, mold release agents, plasticizers, etc.
- Do not use pipet tips that have been sterilized by heating them. The heat may release plasticizer compounds.
- Use Eppendorf LowBind tubes, not siliconized tubes.




Fig. 5. Virus loss in different containers (20mM HEPES, 150 mM NaCl, pH 7.8; initial Ad5 concentration was approximately 6×10^{10} p/mL ≈ 0.02 g/L): (●) LoBind™ protein tube (Eppendorf), (□) Fisherbrand 0.5 mL tube, (■) Glass (National Scientific), (▣) Fisherbrand PCR tube, (○) polycarbonate tube (Eppendorf), (×) prelubricated tube (Costar), (▲) polypropylene tube (ULP), (△) Maxym Recovery (Axygen), (+) silanized glass (National Scientific).

Eppendorf LoBind tubes bind much less viruses then other brand LoBind tubes.

From:
E.I. Trilisky, A.M. Lenhoff: Sorption processes in ion-exchange chromatography of viruses. J. Chromatogr. A 1142 (2007) 2 - 12.

1.4 Recommended procedures

	Protein Identification	Label free relative Quantitation	Silac labeled Relative Quantitation	Dimethyl labeled Relative Quantitation	Absolute Quantitation
FASP	+	+	+	+	+
In Gel Digestion (IGD)	+	+	+		
In Stage Tip (iST)	+	+	+	+	+
On column Dimethyl labeling	-	-	-	+	-
Ovalbumin Standard addition	-	Possible but not advised	-	-	-
Peptide (labeled) Standard addition	-	-	-	-	+
Peptide fractionation	+ (always possible e.g. with High pH RP on μ Columns)				
Measurement	nLC-MS/MS by Proxeon \rightarrow LTQ – Orbitrap XL				
Raw Data handling	MaxQuant (Protein identification and relative quantitation of peptides larger than 7 AA) pNOVO+ (deNovoGUI) for peptide de NOVO sequencing (of peptides smaller than 8 AA)				QualBrowser + MaxQuant
MQ data handling	Perseus (extra filtering of MQ result, statistics and intensity or ratio based clustering of proteins)				
Bioinformatics	<p>Overview of GO tools (but not complete): http://geneontology.org/ First pathway analysis can be done with eg: Reactome (www.reactome.org), PathVisio (http://www.pathvisio.org/), KEGG (http://www.genome.jp/kegg/tool/map_pathway2.html), Biocyc (Subscription needed! Biocyc.org \rightarrow e.g. HumanCyc \rightarrow Change organism and select xxx. Tools \rightarrow Cellular Overview). GO enrichment analysis can be done within: - Internet resources like DAVID (Database for Annotation, Visualization and Integrated Discovery, very easy), or (InterProScan)/PloGO (see below) + WeGO. - Cytoscape (www.cytoscape.org) with plugins like BinGO or ClueGO (or EnrichmentMap). - Use the “R project” environment e.g. PloGO (includes abundance information = more advanced). Interactome studies may benefit from database contained information e.g. in STRING-db.org/ (very easy) but also from Cytoscape plugins like Bionetbuilder or Bisogenet. Alternatives: commercial software like ProteinCenter, Ingenuity or Metacore (these commercial softwares are unfortunately not available at WUR Biochemistry).</p>				

Silac: Stable Isotope Labeling by Amino acids in Cell culture.

FASP: Filter Aided Sample Preparation

MQ: MaxQuant

MaxQuant label free relative quantitation result after analysis with Perseus statistical software. The Figure was prepared in Excel.




Table 1 | Comparison of spin filter, short SDS-PAGE and TFE methods

Method	Protein load	Peptide identifications	Protein identifications
Spin filter	50 μ g	5,369	642
Short SDS-PAGE		4,176	593
TFE		4,663	593
Spin filter	150 ng	86	46
Short SDS-PAGE		298	106
TFE		626	150

Samples of human RKO colon carcinoma cells containing the indicated amounts of protein were prepared in triplicate by the indicated methods and analyzed by reverse phase LC-MS/MS. Peptide identifications are total MS/MS spectrum-to-sequence database matches at 5% false discovery rate; protein identifications are nonredundant identifications with at least two identified peptides and parsimonious protein assembly. Reported values are the means of three technical replicate analyses.

1.5 Sample losses, methods comparison.

Taken from:
Liebler, D. C.; Ham, A. J. L.: Spin filter-based sample preparation for shotgun proteomics. Nature Methods (2009), Volume: 6, Issue: 11, p785-785.

	Yield (%)	Remark
FASP	50	This can be somewhat increased by using multiple digestion steps. (Wisniewski, J.R. and Mann, M. (2012) Consecutive Proteolytic Digestion in an Enzyme Reactor Increases Depth of Proteomic and Phosphoproteomic Analysis. Analytical Chemistry 84(6), 2631-2637)
In-gel digestion (IDG)	20	Purified protein: load a few ug of protein on the gel. Protein extract: more proteins will be identified when you load more protein. Dividing the gel lane into 8 slices means that you should load minimally 16 ug protein on the gel but preferentially (much) more, e.g 50 – 75 ug.
In Stage Tip (iST)	80	

1.6 nLC-MSMS sample necessities.

Hand in minimally 50 ul acidic peptide sample (pH 2 – 4) with a concentration of approximately 10 – 500 ng/ul.

Peptide sample solutions to be measured should not contain detergents (SDS, Tween, Triton, NP etc) or other charged nonvolatile ions like TRIS or phosphate or particles of any kind. (detergent has to be removed by FASP (see section 3.2), by running an SDS gel + doing a In-Gel Digestion or by a specific detergent removal method like with commercial available SDS removal spin columns (e.g. Pierce HiPPR).

Proteins can be transferred into the appropriate (ABC) solution by treatment with a Spin Filter (5 or 10 kDa) filter before doing the Trypsin digestion. Another good and proven method to get rid of impurities is to do FASP (see section 3.2), iST (3.3), an SDS gel-electrophoresis step (see section 4) or at least a C18 uColumn cleanup (see section 1.7). Peptide samples can be desalted/concentrated/transferred by performing the C18 uColumn cleanup as well (see section 1.7).

1.7 General sample cleanup procedures with μ Columns.

Peptide solutions can be concentrated or desalted, and beads (from IP/chromatography/SP3) can be removed using μ Columns (= C18 Stage tip + C18 column material). With this μ Column, there will be no loss of hydrophilic peptides (in contrast to using commercially available micro tips like Zip Tips or StageTips)

1. Prepare your own μ column by:

Use the “cookie cutter” method to cut a small (1.6mm = Gauge 14) piece of a C18 Empore disk (= frit). Transfer the frit to the 200 ul tip with a good fitting plunger and tap it mildly. Add 200 ul methanol to the tip with frit.

Prepare a 50% slurry of LichroprepC18 column material in methanol and add 4 ul of the 50% slurry into the methanol in the tip + frit.

The prepared μ Column can be eluted by hand with a plastic 10 ml syringe, or with the vacuum manifold (connected to a vacuum pump) in lab1. Whatever method you use, do not let the μ Columns run dry.

2. Wash the μ Column ones more with 100 ul MeOH.

3. Equilibrate the μ Column ones with 100 ul 1 ml/l HCOOH in water.

4. Dissolve the sample preferentially in 15 – 200 ul of 1 ml/l HCOOH in water (not containing any AcNi) or in any other aqueous buffer.

5. Add sample(s) to the uColumn(s) (wash gel pieces with an extra 50 or 100 ul 1 ml/l HCOOH in water and add that too) and elute through.

6. Wash the uColumn ones with 100 ul 1 ml/l HCOOH in water.

7. Transfer the uColumn to a new 0.5 ml Eppendorf low binding tube.

8. Manually elute peptides from the C18 Stage tip+ uColumn by adding and eluting with 50 ul 50%AcNi + 50% 1 ml/l HCOOH in water directly into the 0.5 ml Eppendorf low binding tube.

9. The sample is now MALDI-tof ready.

10. For LCMS analysis, reduce the AcNi content by putting the samples in a Concentrator (with open cap) at 45 °C for 2 hours or longer when necessary. The final volume should be below 20 ul. Adjust the sample volume with 1 ml/l HCOOH in water to exactly 50.0 μ l. Sonicate (water bath sonicator) for 5 sec in the hot spot when the sample had been dried completely by accident.

The sample now has <5% AcNi, is pH 3, is particle free and therefore is nLCMS ready. Peptide concentrations can be measured e.g. using Pierce Quantitative Colorimetric Peptide Assay (23275) when the concentration is high enough.

Remark 1: The uColumn cleanup step can also be used when you want to concentrate or combine samples. An alternative option to concentrate or combine peptide samples is the Eppendorf concentrator.

Remark 2: uColumns can be step eluted at pH 10 with 50 ul of a step gradient of AcNi [e.g 5, 8, 11, 18 and 25% AcNi with 20 mM $\text{NH}_4\text{Formate}$ pH10], to **fractionate** peaks according to their high pH affinities for the C18 column. Better: fractionate extensively and mix fractions concatenated: 1+6+11, 2+7+12 etc. (Yang, F., Y. F. Shen, D. G. Camp and R. D. Smith (2012). "High-pH reversed-phase chromatography with fraction concatenation for 2D proteomic analysis." Expert Review of Proteomics 9(2): 129-134)

1.8 About Methionine oxidation

To confirm a peptide sequence containing a Methionine within a known protein, the peptide can be oxidized with 0.1v/v% or 10 μ M H_2O_2 by incubating at 4 °C overnight. This oxidizes every methionine in the sample and shifts the molecular weight up by 16 Da for each methionine.

(Koudelka, T., Dehle, F.C., Musgrave, I.F., Hoffmann, P. and Carver, J.A. (2012) Methionine Oxidation Enhances kappa-Casein Amyloid Fibril Formation. Journal of Agricultural and Food Chemistry 60(16), 4144-4155).

2. Protein determination (BCA).

by the Bicinchoninic Acid (BCA) method (Pierce:#23225): (Sensitivity: 5 – 30 ug protein as measured at 562 nm)

The BCA protein concentration determination method is almost independent on the amino acid composition but cannot be used for samples containing reducing agents (like DTT or TCEP) or chelators (like EDTA).

Alternative: Pierce 660 nm Protein Assay (#22660) with the Ionic detergent Compatibility reagent (#22663) can be used in the presence of maximally 5% SDS and/or 500 mM DTT and is linear between 2.5 and 600 ug BSA or 5 to 1200 ug of Ovalbumin. Compared to BSA, the Ovalbumin response is 54%. This strong protein dependence is the main disadvantage of the Pierce 660 nm Protein Assay.

BCA Solutions:

BCA working solution (BCA ws): 25 ml reagent A + 0.5 ml reagent B
Sample : BCA ws = 1 : 20
Ovalbumin or BSA stock = 1 mg/ml in 50 mM NH_4HCO_3 .

The exact concentration of Ovalbumin or BSA can be directly measured when in 50 mM NH_4HCO_3 at 280 nm:

$$\begin{aligned} C_{\text{Ovalbumin}} &= E_{280} \times 1.30 \quad (\text{mg/ml}) \\ C_{\text{BSA}} &= E_{280} \times 1.50 \quad (\text{mg/ml}) \end{aligned}$$

This is necessary since purified proteins always contain some salts.


Pipet the BCA working solution, water and the Ovalbumin or BSA standard or sample directly in a plastic 1 mL cuvette as described in the table below. Seal tightly with parafilm. Mix and incubate at 21 °C for 60 min or at 37 °C for 30 min. Mix and measure E562.

Example measurement:

Ovalbumin dissolved at 1 mg/ml a absorbance at 280 nm was measured: $E_{280} = 0.618$

Therefore the real concentration was: $0.618 \times 1.30 = 0.803 \text{ mg/ml}$

Volume BCA ws (ul)	Volume Ovalbumin Std (ul)	Volume sample (ul)	Volume H2O (ul)	Calc. amount Ovalbumin (ug)	Measured E562 (1)	Measured E562 (2)
950	0	-	50	0	0.000	0.000
950	10	-	40	8.03	0.186	0.186
950	20	-	30	16.1	0.346	0.351
950	30	-	20	24.1	0.510	0.509
950	40	-	10	32.1	0.653	0.655
950	50	-	0	40.2	0.746	0.730
950	-	20	30			



Linear between 0 – 30 ug:

$$\text{ug protein} = \frac{[E_{562} - 0.0135]}{0.0203}$$

Note: A convenient incubator + absorbance measurement device using 96 well plates is available in lab 4. Decrease all volumes in the table a factor 4 and you can use it to.

3. Gel free protein digestion methods.

When you have a limited amount of protein then you can use a minimal volume (50 μ l) of buffer to do the sonication in and use the iST method (3.1).

3.1 In-StageTip (iST) sample preparation for maximally 5 μ g protein without detergent.

Modified from Ref: Kulak, N. A., G. Pichler, I. Paron, N. Nagaraj and M. Mann (2014). "Minimal, encapsulated proteomic-sample processing applied to copy-number estimation in eukaryotic cells." *Nature Methods* 11(3): 319-324.

Modified (iST) sample prep method:

1. Bring 0.5 mg dried cells or less in a 0.5 ml low binding ep. Wash the cells with 200 μ l 0.1 M Tris pH8 (2*) and suspend them in 100 μ l 8M urea/0.1 M Tris pH 8 or less. Do NOT use SDS or any other detergent. Prepare a homogeneous suspension e.g. by sonication with water bath sonicator (A probe sonicator with a metal tip will contain proteins from all people using it before you adsorbed onto the metal tip).
2. From the cell suspension, pipette 10 μ l protein extract into 200 μ l stage tips with a double C18 membrane and flush through.
3. Add 100 μ l ABC and flush through.
4. + 100 μ l 5% AcNi / ABC and elute through.
5. Reduction: + 20 μ l of 50 mM DTT (7.7 mg/ml ABC). 60°C for 30-60 min, elute through.
6. Alkylation: + 20 μ l of 50 mM AcrylAmide (9 mg/ml 100mM Tris pH8). 21°C in the dark for 30 min (not more), elute through.

Protein cleanup:

7. + 100 μ l ABC, elute through on vacuum station.
8. Do not use the vacuumstation anymore from this point. Put the StageTip into a new 0.5 ml low binding ep. Add 20 μ l Trypsin in ABC (100* dil = 5 ng/ μ l ==> 40 ng) and incubate shaking at rT overnight.
9. Elute through. Add 75 μ l 1 ml/l HCOOH in water to the StageTip, elute through. For high protein loads (> 1 μ g) stage tips may get obstructed here!
10. Add 5 μ l 50% AcNi / 50% 1ml/l HCOOH to the StageTip, elute through, mix low binding ep --> nLCMS ready.

3.2 Filter aided sample preparation (FASP) for maximally 10 μ g protein

Modified from ref: Wisniewski, J. R.; Zougman, A.; Nagaraj, N.; Mann, M.: Universal sample preparation method for proteome analysis. *Nature Methods* (2009) Volume: 6, Issue: 5, Pages: 359-360

1. Sample lysis: Make a concentrated protein sample e.g. by sonicating 20 mg cells in 180 μ l ST-lysis buffer in a 1.5 ml low binding Eppendorf tube.
2. Incubate at 95°C (heating block) for 5 min. Cool down and centrifuge at 13 kRPM for 10 min.
3. Measure the protein content of the supernatant using the BCA method described in section 2. [Protein determination \(BCA\)](#) when necessary. Dilute to 1 μ g/ μ l protein when necessary.
4. Pipet 100 μ g protein sample (e.g. 100 μ l 1 μ g/ μ l) in a new low binding tube, add 0.3 mg of DTT (20 mM) and incubate at 60 C for 1 hour.
5. In a low binding Eppendorf tube (essential!), pipet 450 μ l AA/UT (first, 3.6 mg AA/ml UT) and in the middle add 50 μ l of the reduced protein sample in SDT (50 μ g protein/500 μ l =) and incubate for 1 min while mildly shaking at room temperature.
6. Pipet 10 μ g alkylated protein sample (= 100 μ l) into a Pall 3K omega filter (10-20kDa cutoff, OD003C34) and centrifuge at 13 kRPM for 30 min. ! Do not centrifuge these filters at a higher speed!
7. Add 100 μ l UT solution, gently mix and centrifuge at 13 kRPM for 30 min. Repeat this step 2* to remove all SDS.
8. Add 110 μ l of ABC to the filter unit and centrifuge at 13 kRPM for 30 min.
9. Remove the Pall filter cup from its original micro tube and put it into a clean 2 ml low binding Eppendorf tube.
10. Add 100 μ l ABC containing 1 μ l (0.5 μ g) Trypsin. Centrifuge using the "pulse" button up to 3 kRPM (maximum 3 seconds) and incubate gently shaking at room temperature over night (15 h).
11. Centrifuge the filter units at 13 kRPM for 30 min.
12. For 100 μ l of sample, acidify to pH 2 to 4 by adding 3.5 μ l 10* diluted TFA. The sample is now LCMS ready.

When samples have to be labeled by the Dimethyl labeling method or when you want the absolute maximum amount of peptides in the sample then:

13. Use the acidic ABC eluate for loading onto a C18 μ Column (see section [1.7 General Sample cleanup procedures](#)) and/or Dimethyl labeling (see section [5.1.1 Relative quantitation: on column peptide dimethyl labeling](#)).

3.3 Normal "In solution" trypsin digestion

1. Dissolve 1-10 ug protein in 100 ul 50 mM ABC (pH 8) (10 ug BSA = 0.15 nmol = 5.3 nmol Cys)
2. Add 10 ul 100 mM DTT dissolved into 50 mM ABC (pH 8.5 with NH₃, 1.0 umol). Incubate at 60 °C for 1 hour (max 2 hours).
3. Check the pH. Make pH 8.5 with 10* diluted NH₃. (Usually not necessary!)
4. Add 15 ul 100 mM AcrylAmide dissolved into 50 mM ABC (pH 8.5 with NH₃, 1.5 umol). Incubate at 20 °C in the dark for 0.5 hour exactly.
5. Add 15 ul 125 mM cysteine dissolved into 50 mM ABC to remove the excess Iodoacetamide (1.9 umol) and, for procedure 3.4 only (below), add another 360 ul ABC to decrease the TFE or MeOH concentration to 10%.
6. Add 5 ul trypsin sequencing grade 20* diluted in ABC to 25 ng/ul (125 ng). Incubate gently shaking at room temperature overnight or at 37 °C for 4 h or at 45 °C for 2 h.
7. After digestion, add 2.5 ul 10* diluted TFA to decrease the pH to 2 - 3. Add more 10* diluted TFA when necessary.
8. Perform the sample cleanup with uColumns as described in: 1.7 General sample cleanup procedures with µColumns. In solution digested samples that have not been cleaned by this procedure cannot be measured.

3.4 Methanol and TriFluoroEthanol (TFE) sample preparation method

With Cysteine reduction and CarboxAmidoMethylation (to prevent auto-oxidation) = +57 Da per cysteine.

Ref. Wang, H. X, Qian, W. J.; Mottaz, H. M.;, Smith, R. D.: Development and evaluation of a micro- and nanoscale proteomic sample preparation method. Journal of Proteome Research (2005) 4 (6) p2397-2403

3.5.1. Sonicate in a 2 ml low binding ep 1 mg (or less) of the wet cell sample in 99 ul (or less) of either:

- a. 50 mM ABC (pH 8)
- b. MeOH/50 mM ABC (60:40 v/v)
- c. TFE/50 mM ABC (50:50 v/v)
- d. detergent containing buffer (→ detergent has to be removed by FASP (see section 3.2. Filter aided sample preparation), by running an SDS gel (also see chapter 4 In-Gel Digestion method) or by a specific detergent removal method like with commercial available SDS removal spin columns (e.g. Pierce HiPPR).

3.5.2. Continue with the reduction and alkylation as above under 3.3.2 – 3.3.8.

3.5 Other methods:

SP3: Hughes, C. S., S. Foehr, D. A. Garfield, E. E. Furlong, L. M. Steinmetz and J. Krijgsveld (2014).


"Ultrasensitive proteome analysis using paramagnetic bead technology." Molecular Systems Biology 10(10).

GASP: Fischer, R. and B. M. Kessler (2015). "Gel-aided sample preparation (GASP)-A simplified method for gel-assisted proteomic sample generation from protein extracts and intact cells." Proteomics 15(7): 1224-1229.

Both methods gave poorer results than for iST in my hands but SP3, though rather elaborate, has shown to give good results by others too. E.g.:

Sielaff (2017). "Evaluation Protocols for Proteomic Sample Preparation in the Low Microgram Range." Journal of Proteome Research (in press) with the following figure:

Figure 2




⇒
FASP is OK for ≥ 10 ug protein
iST is OK for ≥ 2 ug protein (when surfactant free)
SP3 works best ≤ 2 ug protein

4. In-Gel Digestion method (IGD)

4.1 General info

Standard protein gel: 12% bisacrylamide (MW 15 – 200 kDa): e.g. Pierce 12% Precise Tris-HEPES precast gels (10 well) (Fisher PX0025202) or from Invitrogen.



Possible internal markers: 2 ug DNA ladder (1 kb plus, Invitrogen 10488-085) added to the sample. Stain after electrophoresis with "Indoine blue" DNA stain (Sigma R325147) as described by:

Guoan Zhang, David Fenyö, and Thomas A. Neubert: Use of DNA Ladders for Reproducible Protein Fractionation by Sodium Dodecyl Sulfate-Polyacrylamide Gel Electrophoresis (SDS-PAGE) for Quantitative Proteomics.


Journal of Proteome Research 7 (2008) 2, p678-686

Suggested external markers to get an indication of the size of a sample protein:

Protein dual color markers (Bio-Rad Precision plus = Cat.# 161-0374) or

Prestained protein MW marker (Fisher PX0026612) 20, 25, 35, 50, 85, 120 kDa

10 ul per well



Comments (by Andrej Shevchenko):

For highest sensitivity, rinse for 60 minutes or more after the gel has been run and fixed. This helps to keep the background transparent during development.

Do not use glutaraldehyde as the sensitizing agent - it is also a protein cross linking agent!

Ref: Shevchenko A., Wilm, M., Vorm, O. and Mann, M. Anal. Chem. T68T, 850-858 (1996). Mass spectrometric sequencing of proteins from silver stained polyacrylamide gels.

4.2 Recommended procedure for CCB or Oriole protein gel staining:

1. With the procedure described on the next page, fixation is not necessary. When you still want to fix the proteins, then do so by incubating the entire gel in: 10% acetic acid / 50% water / 40% methanol. Definitely do NOT use cross-linking agents like glutaraldehyde or formaldehyde!

2. Stain either according to the Colloidal Coomassie Staining method (we use the Colloidal Blue Staining Kit which contains Coomassie G-250 (854.02 g/mol), from Invitrogen ordering# LC6025, 150 euro) described below for visible staining or use the Oriole (Bio-Rad 161-0496 ready for use stain) fluorescent stain to observe bands by UV light.

Colloidal Coomassie Staining can be used to see more than 50 ng of protein, Oriole can go down to the low ng range.

4.2.1 Colloidal Coomassie Preparing Staining Solution

Shake the Stainer B solution before using it. Prepare the solutions fresh (from top to bottom = keep this order and mix after each addition) as described in the table below in a 50 ml (Greiner) tube. Then directly transfer it to the gel in a new square petri dish.

Solution*	ml per gel
Deionized Water	13.75
Methanol	5
Stainer B	1.25
Stainer A	5

*When Stainer A and Stainer B are combined a precipitate may form which will dissolve within 30 seconds.

4.2.2 Colloidal Coomassie Gel staining procedure

1. Shake gel in 25 ml staining solution in a new square petri dish for 1 - 2 hours.

Note: Staining intensity does not vary significantly if left in stain for 3 hours or 12 hours.

2. Decant staining solution and replace with a minimum of 200 ml of deionized water per gel. Shake gel in water for at least 7 hours. The gel will have a clear background after 7 hours in water.

Note: Gels can be left in water for up to 3 days without significant change in band intensity and background clarity.

3. For long-term storage (over 3 days), keep the gel in a 10% Methanol solution at 4°C.

4.3 In Gel Digestion protocol

4.3.1 Remarks

Proteins visible in a SDS gel after Colloidal Coomassie Staining can be measured by Proxeon nLC-LTQ-Orbitrap MS-MSMS provided you read and follow the remarks concerning Keratin in section 1.3 Reduce the amount of Keratins in your samples as much as possible by...

4.3.2 Procedure

1. Sample preparation.

- Make a concentrated protein sample in a suitable buffer e.g. by sonication.
E.g. 50 mg washed cells in 450 μ l SDT-lysis buffer pH 8 (= maximally 5 μ g protein/ μ l) or
2 μ l 2.0 μ g/ μ l BSA in 398 μ l SDT-lysis buffer pH 8 (= 10 ng BSA/ μ l)
- Sonicate and incubate at 95°C (heating block) for 10 min. Cool down to room temperature.
Centrifuge at 13 kRPM for 10 min.
- Pipet 100 μ l of supernatant in an ep and add 25 μ l of Sample loading buffer (e.g. Pierce 39000 or [50% glycerol/50% water + 1 g/l bromophenol blue]). Incubate shaking at 95°C for 10 min. Centrifuge.
The sample will now have circa 1 - 4 μ g protein/ μ l from 50 mg cells or 8 ng BSA/ μ l from 2 μ g BSA
- Apply 2 μ g (one purified protein) to 60 μ g (protein mix) or even **more** protein per sample well.
Each sample well holds maximally 40 μ L (but then you have to pipet very carefully, 30 μ l is more safe).
As a protocol control, use 100 ng of BSA exactly (= 12.5 μ l as prepared above).
Also add to one or two wells: 10 μ l of (pre-stained) protein Molecular weight marker or 2 μ g DNA ladder.

2. Run the gel as described by the manufacturer. For Pierce 12% Precise gels, run at 120V for maximally 55 min = ca 5.5 cm. This is fine for 8 gel slices. If you want to prepare less slices from one sample, then run the gel shorter (ca 5 min per gel slice).

Use gloves. Open the gel cassette with a thin spatula and put it into a large clean square petri dish. Keep the gel covered with the lid as much as possible to prevent extra keratin contamination.

Stain for 1 - 3 hours (see 5.2) and then de-stain with water for 1 to 20 hours. Refresh the water a few times.

Remark: Some (most) in-gel digestion protocols use a 50% acetonitril to wash away the Coomassie but this is not necessary when the sample is measured by LC-MSMS. The Coomassie stain will be released from the protein during the digestion procedure and will elute from the RP column after the peptides at acetonitril percentages higher than 40% but generally is poorly soluble in 1 ml/l HCOOH in water which is used to dissolve the peptides before injection onto the nLC.

3. Cysteines reduction and alkylation.

- Add 25 ml 50 mM NH_4HCO_3 + 35 μ l B-MercaptoEthanol (= 20 mM ME pH 8). Gently shake for 1 h or more at room temperature to reduce all disulfide bridges.
- Wash with water and add 22.5 ml water + 2.5 ml 1M Tris pH 8 + 0.036 g AcrylAmide (= 20 mM AA pH 8).
Incubate at room temperature in the dark while gently shaking for 0.5 hour. Wash with water thoroughly.

4. Gel cutting. (If the gel gets a bit dry and starts jumping around, than add a small drop of water on top of the gel). Cut out the gel bands or slices (1 – 24) and cut them into small pieces of ca. 1 mm². Use a sharp clean scalpel from lab 1 (or a 5 ml white tip or a 1 ml blue tip for higher amounts) on a clean piece of parafilm. Transfer the gel pieces to clean 0.5 ml low binding micro centrifuge tubes. Take a negative (no protein containing gel slice) and positive control (100 ng BSA) as well. A reagent control (containing no gel slice) is not necessary when you use the materials from lab 1 and 2 (rooms 2005 and 2008).

5. Enzymatic digestion.

At this point you may freeze + de-freeze the gel pieces to further increase the Trypsin accessible area.

Add 50 μ l cold freshly prepared Trypsin solution (5 ng/ μ l = 100* diluted into ABC). When there is still some gel piece sticking out of the solution, then add extra ABC (but NO Trypsin) to completely cover the gel pieces.

Preferentially incubate overnight while shaking at room temperature (20 °C) or 4 hours at 37 °C or 2 hours at 45 °C.

6. Extraction of peptides.


- Add 10% TFA up to a pH between 2 and 4 (measure the pH with pH paper, ca. 3.5 μ l is needed per 50 μ l of ABC).
Mix.
- Perform the μ Column cleanup procedure with a C18 μ Column as described in section 1.7 "General sample cleanup procedures with μ Columns". After loading the first peptide extract, wash the remaining gel with 100 μ l 1 ml/l HCOOH in water and add the liquid to the μ Column as well.

The μ Column cleanup step can also be used when you want to concentrate or combine samples. An alternative option to concentrate or combine peptide samples is the Eppendorf concentrator.

5. Quantitation


5.1 Relative quantitation (= Sample compared to Control)

Three ways of relative quantitation can be used. The easiest to do is **label free relative quantitation**. In this case, Samples (to be measured at least in triplo, better 4 times or more) and Controls (the same amounts) are measured separately and relative quantitation takes place after calculation of all peak intensities (by MaxQuant) in each chromatogram. This tricky method became feasible due to the high accuracy and low noise of the Orbitrap's MS measurement and release of the MaxQuant software. The statistical program Perseus can be used to find the really significant differences between samples and controls in the MaxQuant result table. The accuracy strongly depends on the sample preparation reproducibility but generally starts from a factor 10 to find significantly different protein concentrations between sample and control.



Bui, T. P., J. Ritari, S. Boeren, P. de Waard, C. M. Plugge and W. M. de Vos (2015). "Production of butyrate from lysine and the Amadori product fructoselysine by a human gut commensal." *Nat Commun* 6: 10062.


SILAC labeling, essential labeled amino acids like arginin and lysin are added in Light, Intermediate and Heavy form, and these labeled amino acids become incorporated during cell growth. Different samples can therefore be mixed in a very early stage, that is directly after cell lysis and protein determination. Here, protein losses occurring during sample preparation will happen for every labeled protein (light/intermediate and heavy), and therefore will not affect the final result. Silac labeling can be regarded as the most accurate relative quantitation method (Accuracy +/- 30%, you may find significantly different protein concentrations between sample and control from a factor 1.5).



Breast Cancer Cells Exposed to the Phytoestrogen Genistein." *Molecular & Cellular Proteomics* 10(1).

Sometimes Silac labeling is impossible though labeling may be desired. Then, samples can also be labeled on the peptide level. So far, we have obtained best results with **Dimethyl labeling** of N-terminal amines and lysines. This is a reductive alkylation method that uses formaldehyde (CH_2O , CD_2O or $^{13}\text{CD}_2\text{O}$) and cyanoborohydride (NaBH_3CN or NaBD_3CN) and is described in detail below. Disadvantage of labeling at the peptide level is that almost the complete sample preparation has to be done for each sample separately. Different samples are mixed after completion of the sample preparation which may result in a large error. Fortunately, not all proteins will be up- or down regulated due to the stimulus so an internal control should normally be possible (Accuracy: +/- 50%, you may find significantly different protein concentrations between sample and control from a factor 2).


Lu, J., S. Boeren, S. C. de Vries, H. J. F. van Valenberg, J. Vervoort and K. Hettinga (2011). "Filter-aided sample preparation with dimethyl labeling to identify and quantify milk fat globule membrane proteins." *Journal of Proteomics* 75: 340.




5.1.1 Relative quantitation by on column peptide dimethyl labelling protocol

Ref. Boersema, P.J. et al., Nature Protocols (2009) 4, 4, p484-494: Protocol: Dimethyl labeling for relative quantitation.

Reductive amination of NH₂ on N-terminus and lysines with aldehyde:



dΔM = 4 / 8 Da per group



Stock solutions

500 mM NaH ₂ PO ₄ · 1 H ₂ O	(M= 138)	6.9 g/l	3.45 g / 50 ml
500 mM Na ₂ HPO ₄	(M= 142)	7.1 g/l	3.55 g / 50 ml
35% or 20% (vol/vol) formaldehyde in water (CH ₂ O, CD ₂ O or ¹³ CD ₂ O)			
0.6 M cyanoborohydride in water (NaBH ₃ CN or NaBD ₃ CN): (M= 64.8)			3.9 mg / 100 ul

Labeling reagent

Prepare per sample/label:

10 ul 500 mM NaH₂PO₄ · 1 H₂O + 35 ul 500 mM Na₂HPO₄ + 425 ul H₂O (= 50 mM pH 7.5)

+ 2.9 ul 35% or 5 ul of 20% (vol/vol) formaldehyde in water (CH₂O, CD₂O or ¹³CD₂O).

+ 25 ul of 0.6 M cyanoborohydride in water (NaBH₃CN or NaBD₃CN) or 1 mg.

CAUTION Formaldehyde solutions and formaldehyde vapors are toxic, prepare solutions in a fume hood.

CRITICAL Labeling reagent mixtures should be kept at 4 °C and not stored longer than 24 h to ensure labeling efficiency.

Light: CH₂O + NaBH₃CN

Intermediate: CD₂O + NaBH₃CN

Heavy: ¹³CD₂O + NaBD₃CN

On-column stable isotope dimethyl labeling (TIMING: 1 h):

(i) Make C18+ Stage tips by: Use the "cookie cutter" method to put a small (1.6mm) piece of a C18 Empore disk into a plastic 200 ul tip. Do this by cutting the Empore disk with a large metal needle and transfer it to the 200 ul tip with a good fitting plunger from a syringe. Add 200 ul of methanol to this Stage tip. Add 5 ul of a 50% LichroprepC18 slurry in methanol to the Stage tip to create a C18+ Stage tip.

(ii) Wash the C18+ Stage tip column with 200 ul of methanol.

(iii) Condition the C18+ Stage tip column with 100 ul of 1 ml/l HCOOH in water.

(iv) Load the acidified peptide sample in water (maximally 10 ug peptide, detergent free) on a C18+ Stage tip column.

(v) Wash the C18+ Stage tip column with 100 ul of 1 ml/l HCOOH in water.

(vi) Put 2 ml Eppendorf tubes under the columns to catch reagent that runs through the column.

(vii) In 10 min. time (not faster), flush each of the C18+ Stage tip columns with 100 ul of the respective labeling reagent (light, intermediate or heavy).

CRITICAL STEP To allow for complete labeling, make sure that Step vii takes at least 10 min.

(viii) Wash the C18+ Stage tip columns with 200 ul of 1 ml/l HCOOH in water.

(ix) Add 10 ul 1 M Tris to the eps under the columns to destroy remaining reagent.

Peptide recovery

Manually (= with a syringe) elute and collect the labeled samples in new 0.5 ml low binding eps from the C18 Stage tip columns with 50 ul of 50% AcNi/50% 1 ml/l HCOOH in water.


CRITICAL STEP When performing the protocol for the first time or with a new sample, it is advised to check the labeling efficiency and sample amounts by measuring a fraction of the sample by LC-MS before mixing differentially labeled samples.

For LCMS analysis, reduce the AcNi content by putting the samples in a Concentrator (with open cap) at 45 °C for 2 hours or longer when necessary. The final volume should be below 20 ul. Adjust the sample volume with 1 ml/l HCOOH in water to exactly 100.0 ul. Sonicate (water bath sonicator) for 5 sec in the hot spot when the sample had been dried completely by accident.

The sample now has <5% AcNi, is pH 3 and is particle free. Mix the samples in the appropriate ratio (light/intermediate/heavy) to make them nLCMS ready.

5.2 Absolute quantitation


- a.** Absolute quantitation on a limited number of proteins can be done by making **calibration curves with synthetic peptides that preferentially contain 1 ¹³C labeled amino acid**, e.g. the C-terminal K or R in case of a tryptic peptide. To do this, the (HPLC or NMR quantified) labeled peptide has to be added to the sample to generate the calibration curve under exactly the same conditions as the sample measurements. Later, the labeled peptide can also serve as an internal standard. This way of absolute quantitation generally gives a good accuracy (+/- 30%).



- b.** Less stringent but more practical when more proteins have to be quantitated, is the method of quantitation with respect to a single added internal standard like BSA (e.g. Pierce **PI23210** 2.0 mg/ml) or Ovalbumin. This rough absolute quantitation “**relative to an internal standard**” is accurate within a factor of 4.

BSA can be replaced by Waters MassPREP Digestion Standard Mix 1 [186002865] (a 4 protein mix, €200/set) or Sigma UPS2, a quantified set of 48 human proteins at different concentrations (€900/set).

Lahtvee, P.J.: Absolute quantification of protein and mRNA abundances demonstrate variability in gene-specific translation efficiency in yeast. *Cell Syst.* 4, 495–504.e5 (2017).





- c.** Somewhat less accurate but very easy to implement is the “**intensity based absolute quantification (iBAQ)**” method (Schwanhauser et al. *Nature* 2011, 473, 7347, P337-342) that uses the Total peak intensity as determined by MaxQuant for each protein and corrects that for the number of measurable peptides (=number of tryptic peptides of 7-30 amino acids long without missed cleavages).

The error in the iBAQ method is between a factor 2 (as shown to the left) and 10 as shown in Log10 iBAQ vs Log [Protein] plot.

The Log-Log figure comes from a presentation by Selbach during the MaxQuant summerschool 2011. It has been published as Fig. S8 in Schwanhauser 2011 mentioned above.

Absolute protein levels: intensity-based absolute quantification (iBAQ)




- d.** For very large datasets obtained by a Q-exactive or Fusion (>12.000 peptides), the absolute amount of a protein per cell can be calculated by relating the MS signal intensity from a protein to the total MS signal intensities of all histones observed. Error = a factor 2 only (?). Wisniewski (2014): A ‘proteomic ruler’ for protein copy number and concentration estimation without spike-in standards. MCP

- e.** Addition of labeled peptides produced by either Cell free expression systems (Takemori: MEERCAT Multiplexed Efficient Cell Free Expression of Recombinant QconCATs For Large Scale Absolute Proteome Quantification. *Mol Cell Proteomics.* 2017, 16(12):2169-2183. SB ca €5000) or Ecoli (Beynon, R. J., Doherty, M. K., Pratt, J. M., and Gaskell, S. J. (2005) Multiplexed absolute quantification in proteomics using artificial QCAT proteins of concatenated signature peptides. *Nat. Methods* 2, 587–589 SB ca €3000).

6. Phosphopeptide enrichment methods (S and T only).


Serine, Threonine: TiO₂ will yield more mono phosphorylated S + T peptides than IMAC.
IMAC will yield relatively more doubly phosphorylated S + T peptides than TiO₂
Tyrosine: Use specific anti-bodies (e.g. Cell Signaling PhosphoScan kit P-Tyr-1000 #8803S)
RP LC coupled to electrospray ionization will mainly show singly and doubly phosphorylated peptides.
Maldi ionization may give more multiply phosphorylated peptides but also Maldi is less sensitive to multiple phosphorylated peptides.



B. Bodenmiller, L. N. Mueller, M. Mueller, B. Domon and R. Aebersold Reproducible isolation of distinct, overlapping segments of the phosphoproteome. *Nature Methods* (2007) 4, 3, p231-237.



The SCX/IMAC Enrichment Approach for Phosphoproteomics



A protocol often used:

B. Macek, M. Mann and J. V. Olsen: Global and Site-Specific Quantitative Phosphoproteomics: Principles and applications. *Annual Review of Pharmacology and Toxicology* (2009) 49, p199-221.

Sample: 7.5x10⁷ human K562 human chronic myelogenous leukemia cells, 4mg lysate
Protocol: Villen, J, and Gygi, SP, *Nat Prot*, 2208, 3, 1630-1638.
Lysis: 8M urea, 75mM NaCl, 50 mM Tris pH 8.2, phosphatase inhibitors
SCX: PolyLC - Polysulfoethyl A 9.4 mm X 200mm, elute: 0-105mM KCl, 30% Acn.
IMAC: Sigma - PhosSelect Fe IMAC beads, bind: 40% Acn, 0.1% formic acid, elute: 500 mM K₂HPO₄ pH 7
MS/MS: Thermo Fisher Orbitrap XL, high-res MS1 scans in the Orbitrap (60k), Top-8 fragmented in LTQ, exclude +1 and precursors w/ unassigned charges, 20s exclusion time, precursor mass error +/- 10 ppm

6.1 Phosphopeptides (S, T) enrichment by Titanium Dioxide (TiO₂) Chromatography.

6.1.1. Enrichment with normal TiO₂ in uColumns.

Aryal, U. K. and A. R. S. Ross (2010). "Enrichment and analysis of phosphopeptides under different experimental conditions using titanium dioxide affinity chromatography and mass spectrometry." *Rapid Communications in Mass Spectrometry* 24(2): 219-231.

Zhou, H. J., T. Y. Low, M. L. Hennrich, H. van der Toorn, T. Schwend, H. F. Zou, S. Mohammed and A. J. R. Heck (2011). "Enhancing the Identification of Phosphopeptides from Putative Basophilic Kinase Substrates Using Ti (IV) Based IMAC Enrichment." *Molecular & Cellular Proteomics* 10(10).

Humphrey, S. J., S. B. Azimifar and M. Mann (2015). "High-throughput phosphoproteomics reveals in vivo insulin signaling dynamics." *Nat Biotechnol* 33(9): 990-995.

Modified procedure:

1. To 90 μ L of the centrifuged digested peptide solution (Mix after each addition):
2. Add 10 μ L TriFluoroEthanol
3. Add: 114 μ L Acetonitril → samples may become slightly turbid
4. Add: 14 μ L TriFluoroAcetic acid → 50% ACN, 6% TFA
5. Mix peptide solutions at room temperature for 1 min in a ThermoMixer at 300 rpm.
6. Centrifuge at high speed ($\geq 16,000$ xg for 15 min).
7. Prepare ca 5 mm uColumns with 10 μ L TiO₂ beads slurry (50% beads in 100% MeOH, GL Sciences #5010-21315) pipetted into 200 μ L AcNi using a 200 μ L tip fitted with a C8 filter. With 5 mg TiO₂ beads you can use 25 to 50 μ g of peptide solution to be enriched.
8. Wash the TiO₂ column with 100 μ L ultrapure AcNi (e.g. HPLC gradient grade).
9. Equilibrate the TiO₂ column with 200 μ L Loading buffer (80% ACN, 6% TFA).
10. Add sample to the TiO₂ column and slowly elute in 5 min at 18 bar Hg.
11. Non-specifically bound peptides are washed from the TiO₂ beads with:
12. 1* 200 μ L Loading buffer (80% ACN, 6% TFA) and
13. 2* 200 μ L Wash buffer (60% ACN, 1% TFA). In the **last step, elute until the column runs just dry.**
14. Elute Phosphopeptides into low-binding tubes slowly with 50 μ L Elution buffer (40% ACN, 15% NH₄OH (25%, HPLC grade) added, prepare immediately before use). Pulsed (20 μ L, 30 μ L) with 5' on shaker in between pulses.
15. Concentrate samples in a concentrator for 30 min at 45°C.
16. Add 10 μ L 10% TFA or more until the samples are acidic.
17. Perform the general C18 uColumn cleanup as in 1.7.

6.1.2. Enrichment with magnetic Ti⁴⁺ beads


Vu, L. D., E. Stes, M. Van Bel, H. Nelissen, D. Maddelein, D. Inze, F. Coppens, L. Martens, K. Gevaert and I. De Smet (2016). "Up-to-Date Workflow for Plant (Phospho)proteomics Identifies Differential Drought-Responsive Phosphorylation Events in Maize Leaves." *J Proteome Res.*

Alternatively use a magnetic resin for more convenient phosphopeptide enrichment. Here a MagReSyn Ti-IMAC microspheres (ReSyn Biosciences) is used which is **stable for 6 months** when kept cold. Beads can be obtained from Mark Roosjen who tested this method.

1. Mix the beads thoroughly and pipet 50 μ L = 1 mg of MagReSyn Ti-IMAC microspheres (ReSyn Biosciences) into a 2 ml low binding ep. Clear the suspension by putting it on the magnet and remove the supernatant.
2. Wash the beads with 200 μ L of 70% ethanol, and remove the sup.
3. Just before loading with sample, wash the beads with 100 μ L of loading solvent (80% acetonitrile, 5% TFA). Remove the sup.
4. Dissolve the desalted and dried trypsin phosphopeptides containing digest into 100 μ L of loading solvent (80% acetonitrile, 5% TFA).
5. Incubated with 1 mg of MagReSyn Ti-IMAC washed and equilibrated beads for 20 min at room temperature.
6. Wash once with 100 μ L wash solvent 1 (80% acetonitrile, 1% TFA, 200 mM NaCl).
7. Wash twice with solvent 2 (80% acetonitrile, 1% TFA).
8. Elute the bound phosphopeptides with 80 μ L of a fresh 1% NH₄OH solution.
9. Immediately acidify the eluate to pH ca 3 with 10% TFA. Check the pH by putting 0.1 μ L on a piece of pH paper.

6.2 Phosphopeptide (S,T) sample preparation by sequential elution from IMAC.

Uses the Pierce Fe-NTA Phosphopeptide Enrichment Kit: PX0088300



Ref: Thingholm, T.E., Jensen, O.N., Robinson, P.J. and Larsen, M.R. (2008) SIMAC (sequential elution from IMAC), a phosphoproteomics strategy for the rapid separation of mono phosphorylated from multiply phosphorylated peptides. *Molecular & Cellular Proteomics* 7(4), 661-671.

This protocol uses a step elution procedure to separate multi-phosphorylated peptides from mono- and di-phosphorylated peptides. Mono- and di-phosphorylated peptides are made free of non-phosphorylated peptides by using TiO₂ columns.

A. Sample Preparation

1. Perform a tryptic digestion, preferentially using the FASP protocol (see [3.2. Filter aided sample preparation \(FASP\)](#)).
2. Quick dry the peptide samples in a rotary vacuum concentrator.
3. Resuspend the dried sample to a concentration of 1-10 ug/ul in Loading Buffer (= 0.1% TFA, 50% AcNi).
4. Wash the Fe-NTA spin column (which contains 200 ul slurry) with 2 * 500 ul of Loading Buffer.
5. Add maximally 5 mg peptides to a washed Fe-NTA spin column and incubate for 20 minutes at room temperature with end-over-end rotation. Remove bottom tab of the column. Place column in a microcentrifuge tube.
6. Centrifuge column at 2,000 rpm for 2 minutes. Collect the flow-through for analysis = **non-bound**. Transfer the column to a new collection tube.

B. Wash and elute from IMAC

Wash with 150 ul Loading buffer = 0.1% TFA, 50% AcNi
 Wash 2* with 150 ul Wash buffer = 1 % TFA, 20% AcNi
 Wash with water.

Flow through = Sample **A**
 Flow through = Sample **B**

C. Elute with 100 ul ammonia water, pH 11.3 (10 ul of 25% ammonia solution in 490 ml of MilliQ water) and decrease the pH to around 3 with 10% TFA. Flow through = Sample **C**

Samples non-bound, A and B contain non-phosphorylated peptides and abundant mono-phosphorylated peptides. Sample C contains by far most mono and multi-phosphorylated peptides.

Dry sample C in a vacuum evaporator at 45 °C for 2 - 4 h. Sample C can be dissolved into 100 ul 1 ml/l HCOOH to be analyzed by LCMS.

When you want to see everything: Titanium Dioxide (TiO₂) Chromatography for Samples unbound, A and B.

Dry all samples in vacuum evaporator at 45 °C for 2 - 4 h. Unbound, A and B.

Perform the TiO₂ chromatography step to samples A and B obtained at step 7.2C to get rid of non-phosphorylated peptides as described under 7.1 Phosphopeptides enrichment by Titanium Dioxide (TiO₂) Chromatography.

Dissolve dried samples (unbound, A and B) in 100 (A, B) or 1000 (unbound) ul 1 ml/l HCOOH. Samples are now nLCMS ready.