

[protocol] BCA method

woensdag 13 mei 2015 11:05

Kit is in sample fridge

- Make sure the samples are transparent (filter if needed with syringe filter 0.2 μ m)
Estimate the calibration curve range: e.g. 0.1% max for Secondary emulsions nonadsorbed? = 1 g/L max

- Make sure you have calibration samples ready (fresh or defrozen):

Stock: 10 g/L Citric acid buffer (10 mM; pH 3)

Conc (g/L)	0	0.2	0.4	0.6	0.8	1
Vstock (uL)	0	20	40	60	80	100
Vbuffer (uL)	1000	980	960	940	920	900
Vtotal	1 mL	1 mL	1 mL	1 mL	1 mL	1 mL




Or Stock: 1.0 ml of a standard solution consisting of 1.0 g/L BSA in 0.15 M NaCl with 0.05% sodium azide as a preservative.. Lower volume calibration curve:


Conc (g/L)	0	0.2	0.4	0.6	0.8	1
Vstock (uL)	0	0.2*300=60	0.4*300=120	0.6*300=180	0.8*300=240	300
Vbuffer (uL)	400	340	280	220	160	0
Vtotal	0.2 mL					

60+120+180+240+300=900 mL

DEZE EPJES INVRIEZEN!!!

- Prepare fresh bca working reagent:
50 parts of A plus 1 part B, shelf life of a day
so for 12*3+6=42 tubes you need for example 50 ml to be sure: 50ml A plus 45/50=1 ml B
- Take 50ul from the calibration solutions (0 tot 1 g/L calibration curve; 50 uL sample; als je verwacht dat je sample hogere concentratie heeft moet je hem vooraf verdunnen! Niet minder toevoegen)
Can be frozen in between.
- Dilute "samples" 5x and 10x in buffer: 100 uL in 400 / 900 uL buffer respectievelijk, vortex!
- Take twice 50ul of the undiluted and diluted samples (A, B, C) and add them to another eppendorf tube
- Add 1ml working reagent to all these epjes (A, B, C, calibration curves).
- Incubate 2h at room temperature, or 30 min at 37°C in shaking device
- Cool to fridge 5min, then 10 min room
- Measure absorbance at Fixed wavelength at 562 nm: in 1ml cuvettes
 - Blanc is the zero concentration, measure this first
 - Start timer and measure calibration curve
 - Measure samples
 - Measure calibration curve
 - Write down total time of measuring: min (should be <10 min!)
- Turn off spectrophotometer
- Dispose the reagent in the sink with a lot of rinsing

New:	Old:
 bca1bul	 BCA method
 BCA checklist	



Concentration (g/L)	absorbance at 562 nm
0	
0.2	
0.4	
0.6	
0.8	
1	

Primary emulsion	Abs 562 nm	duplo
Sample before (pH3) A		
Sample before (pH3) B		
Sample before (pH3) C		
Sample after (pH7) A		
Sample after (pH7) B		
Sample after (pH7) C		

Secondary	Abs 562 nm	duplo
Sample before (pH3) A		
Sample before (pH3) B		
Sample before (pH3) C		
Sample after (pH7) A		
Sample after (pH7) B		
Sample after (pH7) C		

Tertiary emulsion	Abs 562 nm	duplo
Sample before (pH3) A		
Sample before (pH3) B		
Sample before (pH3) C		

Sample after (pH7) A		
Sample after (pH7) B		
Sample after (pH7) C		

Concentration (g/L)	absorbance at <u>562 nm</u> Duplo at the end
0	
0.2	
0.4	
0.6	
0.8	
1	