

# Catalogue 2024

---

**primediagnos**tics

Droevendaalsesteeg 1  
P.O. Box 16  
6700 AA Wageningen  
The Netherlands

phone : +31 317 480 613  
e-mail : [primediagnos@wur.nl](mailto:primediagnos@wur.nl)  
url : [www.primediagnos.com](http://www.primediagnos.com)  
webshop: [shop.wur.nl/primediagnos](http://shop.wur.nl/primediagnos)



## Table of Contents

Antisera Product information .....	1
General information on antisera .....	1
Storage information and stability .....	1
Information concerning the quality of the antisera .....	1
Antisera to plant-pathogenic bacteria .....	1
Virus inocula .....	1
Orders and invoicing .....	2
Contact persons for information .....	2
Abbreviations .....	2
Antisera Virus .....	3
Antisera Bacteria and Fungi .....	4
Prime Diagnostics Terms of Delivery .....	5
Workflow overview .....	6
Schematic overview of the ELISA procedure .....	6
Schematic overview of the Luminex procedure .....	7
Schematic overview of the TaqMan procedure .....	8
ELISA and Luminex xMAP comparison .....	9
Appendix 1. Accessing webshop .....	10

## General information on antisera

- Antisera that are not on the list can be produced on request. Please inquire.
- Custom-made conjugates with fluorescent labels or AP can be produced upon request.
- Protocols are available at [www.primediagnosics.com](http://www.primediagnosics.com).
- The AP-conjugate and coating should be diluted 1000x.
- The coating and alkaline phosphatase (AP) conjugates are suitable for use in DAS-ELISA only.

## Storage information and stability

- After long time storage at -20 °C avoid repeated freezing and thawing, as this can influence the product quality considerably.
- Coatings and AP-conjugates should preferably be stored at 4 °C. At 4 °C antisera are stable for at least two years.
- Controls and IgG should be stored in aliquots at -20 °C. Avoid repeated freezing and thawing. When treated properly these reagents will remain stable for at least two years.
- For long time storage (periods longer than one year) store at -20 °C, prepare aliquots if necessary.
- Lateral flow devices, packed in their original unopened (sealed) envelopes can be stored at room temperature. They will remain stable until expiration date. Lateral flow devices, packed in their original opened envelopes will remain stable for two weeks.
- Luminex reagents should be stored at 4°C. They will remain stable for two years.
- Lyophilized positive control (VPC) contain lyophilized leaf material. The lyophilized VPC remain stable for 1 year when stored at 4 °C. A prepared VPC can be stored in aliquots of 0.2 mL and is stable for 1 month at -20 °C and at 4 °C for 1 day.

## Information concerning the quality of the antisera

- Antisera are biological products and differences in reactivity between batches are apparent. Prime Diagnostics exercises strict quality control and guarantees the reactivity as indicated on the batch quality sheet (if applicable).
- Negative controls are optimized for use in different crops and are not suitable for determination of the selectivity of the assay.
- Positive controls are qualitative and cannot be used for quantification!
- The criteria for reactivity of ELISA reagents for bacterial detection are an OD<sub>405</sub> reading of the 10 times diluted positive control (end concentration of 10<sup>7</sup> cells/mL) of OD<sub>405</sub> > 1.0 after 30 minutes of substrate incubation.
- The criterion for reactivity of ELISA reagents for virus detection is an OD<sub>405</sub> reading of the 10 times diluted positive control of at least 0.5 OD after 30 minutes of substrate incubation. The value of the OD<sub>405</sub> may however vary depending on the ELISA reader used.
- Within our quality control procedure, our antisera are tested for lack of reactivity against a variety of pathogens and/or plant compounds known to occur in combination with the given pathogen. However, within the context of any given test, the possibility of false-positive results with unrelated organisms or plant or matrix substances should always be considered.
- Within our quality control procedure, our antisera are tested for specificity against a wide variety of currently known and recognized pathogenic isolates, strains and patho-types of the subsequent pathogens. However, new pathogenic or non-pathogenic strains and isolates of a given pathogen may emerge over time. Prime Diagnostics cannot guarantee that such a new strain or isolate will be detected with equal efficiency.

## Antisera to plant-pathogenic bacteria

- Antisera labeled with different conjugates are available on request. Prices can be supplied upon request.
- For all bacteria non-infectious positive controls are available for use in DAS-ELISA or IIF. These controls have a concentration of 10<sup>8</sup> cells/mL.
- Negative controls consisting of non-infectious bacteria or plant extract are available on request.
- The dilution of the IgG for IIF is dependent on the assay used and laboratory circumstances. It is recommended to determine the working dilution of the reagents under on-site conditions.
- The IgG is suitable for use in indirect immunofluorescence (IIF) detection.

## Virus inocula

- Freeze-dried virus particles, dried plant material (to be reconstituted in buffer) or purified virus for inoculating purposes can be obtained on request. The amounts are generally sufficient for small-scale screening experiments.

## Orders and invoicing

phone : +31 317 480 613  
 e-mail : [primediagnosics@wur.nl](mailto:primediagnosics@wur.nl)  
 url : [www.primediagnosics.com](http://www.primediagnosics.com)  
 webshop : [shop.wur.nl/primediagnosics](http://shop.wur.nl/primediagnosics)

## Contact persons for information

Sales	<a href="mailto:primediagnosics@wur.nl">primediagnosics@wur.nl</a>		
General information	José van Beckhoven	(+31) 317 480 603	<a href="mailto:jose.vanbeckhoven@wur.nl">jose.vanbeckhoven@wur.nl</a>
ELISA	José van Beckhoven	(+31) 317 480 603	<a href="mailto:jose.vanbeckhoven@wur.nl">jose.vanbeckhoven@wur.nl</a>
	Anke Clerkx	(+31) 317 480 663	<a href="mailto:anke.clerkx@wur.nl">anke.clerkx@wur.nl</a>
IIF	Margarit de Klein	(+31) 317 480 669	<a href="mailto:margarit.deklein@wur.nl">margarit.deklein@wur.nl</a>
Inocula	Petra van Bekkum	(+31) 317 486 615	<a href="mailto:petra.vanbekkum@wur.nl">petra.vanbekkum@wur.nl</a>
LAMP	Stefan Aanstoot	(+31) 317 486 615	<a href="mailto:stefan.aanstoot@wur.nl">stefan.aanstoot@wur.nl</a>
	Yvonne Griekspoor	(+31) 317 480 763	<a href="mailto:yvonne.griekspoor@wur.nl">yvonne.griekspoor@wur.nl</a>
	Willem Stol	(+31) 317 480 561	<a href="mailto:willem.stol@wur.nl">willem.stol@wur.nl</a>
Lateral Flow Devices	Jan Bergervoet	(+31) 317 480 836	<a href="mailto:jan.bergervoet@wur.nl">jan.bergervoet@wur.nl</a>
	Antje de Bruin	(+31) 317 480 631	<a href="mailto:antje.debruin@wur.nl">antje.debruin@wur.nl</a>
Luminex	Jan Bergervoet	(+31) 317 480 836	<a href="mailto:jan.bergervoet@wur.nl">jan.bergervoet@wur.nl</a>
	Sander Grapendaal	(+31) 317 484 462	<a href="mailto:sander.grapendaal@wur.nl">sander.grapendaal@wur.nl</a>
	Henry van Raaij	(+31) 317 480 943	<a href="mailto:henry.vanraaij@wur.nl">henry.vanraaij@wur.nl</a>
Marketing and sales	Willem Stol	(+31) 317 480 561	<a href="mailto:willem.stol@wur.nl">willem.stol@wur.nl</a>
	Patrick de Leeuw	(+31) 317 482 910	<a href="mailto:patrick.deleeuw@wur.nl">patrick.deleeuw@wur.nl</a>
Packaging and shipping	Fianne Lubbers	(+31) 317 486 614	<a href="mailto:fianne.lubbers@wur.nl">fianne.lubbers@wur.nl</a>
	Antje de Bruin	(+31) 317 480 631	<a href="mailto:antje.debruin@wur.nl">antje.debruin@wur.nl</a>
Prime Support	José van Beckhoven	(+31) 317 480 603	<a href="mailto:jose.vanbeckhoven@wur.nl">jose.vanbeckhoven@wur.nl</a>
	Jan Bergervoet	(+31) 317 480 836	<a href="mailto:jan.bergervoet@wur.nl">jan.bergervoet@wur.nl</a>
TaqMan	Sander Grapendaal	(+31) 317 484 462	<a href="mailto:sander.grapendaal@wur.nl">sander.grapendaal@wur.nl</a>
	Henry van Raaij	(+31) 317 480 943	<a href="mailto:henry.vanraaij@wur.nl">henry.vanraaij@wur.nl</a>
Invoicing	Marieke Forch	(+31) 317 480 640	<a href="mailto:marieke.forch@wur.nl">marieke.forch@wur.nl</a>

## Abbreviations

IIF	indirect immunofluorescence	CO	coating
DAS	double antibody sandwich	OD <sub>405</sub>	optical density at 405 nm
ELISA	enzyme linked immuno sorbent assay	mL	milliliter
LFD	later flow device	PC	positive control
AP	alkaline phosphatase	pv.	pathovar
VPC	lyophilized positive control	subsp	subspecies

The antisera listed below are available for ELISA and Luminex format, LFD format on request. Ordering information and details of these formats can be found in the webshop of Prime Diagnostics: [shop.wur.nl/primediagnostics](http://shop.wur.nl/primediagnostics).

Pathogen (virus)	Acronym	Pathogen (virus)	Acronym
Alfalfa mosaic virus	AMV	Melon necrotic spot virus	MNSV
Alstroemeria carlavirus	AICV	Mirafiori lettuce big vein virus	MiLBVV
Alstroemeria flower-banding virus	AIFBV	Odontoglossum ringspot virus	ORSV
Alstroemeria mosaic virus	AIMV	Onion yellow dwarf virus	OYDV
Alternanthera mosaic virus	AltMV	Pea early-browning virus	PEBV
Andean potato latent virus/Andean potato mild mosaic virus	APLV/ APMMV	Pea seed-borne mosaic virus	PSbMV
Andean potato mottle virus	APMoV	Pelargonium flower-break virus	PFBV
Apple chlorotic leafspot virus	ACLSV	Pelargonium line pattern virus	PLPV
Apple stem grooving virus	ASGV	Pepino mosaic virus	PepMV
Arabis mosaic virus	ArMV	Pepper mild mottle virus	PMMoV
Bean common mosaic necrosis virus	BCMNV	Plum pox virus ('sharka')	PPV
Bean common mosaic virus	BCMV	Potato leafroll virus	PLRV
Bean yellow mosaic virus	BYMV	Potato virus A	PVA
Beet western yellows virus	BWYV	Potato virus M	PVM
Bell pepper mottle virus	BePMV	Potato virus S	PVS
Calibrachoa mottle virus	CbMV	Potato virus V	PVV
Carnation etched ring virus	CERV	Potato virus X	PVX
Carnation latent virus	CLV	Potato virus Y	PVY
Carnation mottle virus	CarMV	Prune dwarf virus	PDV
Carnation necrotic fleck virus	CNFV	Shallot yellow stripe virus	SYSV
Carnation ringspot virus	CRSV	Squash mosaic virus	SqMV
Carnation vein mottle virus	CVMV	Strawberry latent ringspot virus	SLRSV
Cherry leafroll virus	CLRV	Streptocarpus flower-break virus	SFBV
Chrysanthemum virus B	CVB	Tobacco mild green mosaic virus	TMGMV
Cucumber green mottle mosaic virus	CGMMV	Tobacco mosaic virus	TMV
Cucumber mosaic virus	CMV	Tobacco ringspot virus	TRSV
Cymbidium mosaic virus	CymMV	Tomato aspermy virus	TAV
Dasheen mosaic virus	DsMV	Tomato black ring virus	TBRV
Hosta virus X	HVX	Tomato brown rugose fruit virus	ToBRFV
Impatiens necrotic spot virus	INSV	Tomato bushy stunt virus	TBSV
Iris yellow spot virus	IYSV	Tomato mosaic virus	ToMV
Kalanchoe mosaic virus	KMV	Tomato ringspot virus	ToRSV
Leek yellow stripe virus	LYSV	Tomato spotted wilt virus	TSWV
Lettuce big vein associated virus	LBVaV	Tulip Virus X	TVX
Lettuce mosaic virus	LMV	Zucchini yellow mosaic virus	ZYMV

The antisera listed below are available for IIF, other formats on request. Ordering information and details can be found in the webshop of Prime Diagnostics: [shop.wur.nl/primediagnosics](http://shop.wur.nl/primediagnosics)

Pathogen (bacteria)	Acronym
<i>Curtobacterium flaccumfaciens</i> pv. <i>flaccumfaciens</i>	Cff
<i>Curtobacterium flaccumfaciens</i> pv. <i>oortii</i>	Cfo
<i>Clavibacter insidiosus</i>	Ci
<i>Clavibacter michiganensis</i> subsp. <i>michiganensis</i>	Cmm
<i>Clavibacter michiganensis</i> subsp. <i>sepedonicus</i>	Cms
<i>Dickeya chrysanthemi</i> ( <i>Erwinia chrysanthemi</i> )	Ech
<i>Erwinia amylovora</i>	Eam
<i>Pectobacterium atrosepticum</i>	Eca
<i>Pseudomonas chichorii</i>	Pc
<i>Pseudomonas syringae</i> pv. <i>lachrymans</i>	Psi
<i>Pseudomonas syringae</i> pv. <i>mors-prunorum</i>	Psm
<i>Pseudomonas savastanoi</i> pv. <i>phaseolicola</i>	Psph
<i>Pseudomonas syringae</i> pv. <i>pisi</i>	Pspi
<i>Pseudomonas syringae</i> pv. <i>porri</i>	Pspo
<i>Pseudomonas savastanoi</i> pv. <i>savastanoi</i>	Pssa
<i>Pseudomonas syringae</i> pv. <i>syringae</i> (strain specific)	Pssy
<i>Pseudomonas syringae</i> pv. <i>tomato</i>	Pst
<i>Rhodococcus fascians</i>	Rhf
<i>Ralstonia solanacearum</i>	Rsol
<i>Xanthomonas axonopodis</i> pv. <i>dieffenbachiae</i>	Xcd
<i>Xanthomonas phaseoli</i> pv. <i>phaseoli</i>	Xcph
<i>Xanthomonas citri</i> pv. <i>fuscans</i>	Xcphf
● <i>Xanthomonas axonopodis</i> pv. <i>begoniae</i>	Xcb
<i>Xanthomonas campestris</i>	Xccam
<i>Xanthomonas axonopodis</i> pv. <i>vesicatoria</i>	Xcv
<i>Xanthomonas axonopodis</i> pv. <i>fragariae</i>	Xf
<i>Xanthomonas hyacinthi</i>	Xch
<i>Xanthomonas hortorum</i> pv. <i>pelargonii</i>	Xcp
<i>Xanthomonas arboricola</i> pv. <i>pruni</i>	Xpru

The antiserum listed below is available for LFD, other formats on request. Ordering information and details can be found in the webshop of Prime Diagnostics: [shop.wur.nl/primediagnosics](http://shop.wur.nl/primediagnosics)


Pathogen (fungi)	Acronym
<i>Phytophthora</i>	Phyto

### Special terms and conditions of Wageningen Plant Research (WPR) concerning Prime Diagnostics Products

- These “Special terms of delivery” of Wageningen Plant Research (WPR) and the “General Terms and Conditions” of Wageningen University & Research will apply to all offers and all agreements between WPR and client concerning Prime Diagnostics Products. In case of conflict between both terms, these “Special terms and conditions” takes precedence over the “General Terms and Conditions”.
- Only webshop and e-mailed orders that have been accepted by WPR will be handled. Orders by phone are not handled.
- Complaints regarding deficiencies of the products delivered by WPR should be deposited within 2 months after the purchase and should be supported by relevant test results obtained with the standard protocol of Prime Diagnostics. In the event that WPR declares a complaint to be founded, WPR will be exclusively obliged to effect performance as agreed upon as yet or to refund the purchase price paid, at WPR's exclusive discretion.
- All agreements between WPR and the clients are governed by Dutch Law only. All disputes shall be handled exclusively by the competent court in The Hague, The Netherlands.
- All products are supplied under the condition that they are for the exclusive use by the client. They may not be sold, integrated into other commercial applications nor handed over to third parties without distribution or supply and license agreement.
- Antisera are biological products and differences in reactivity between batches may occur. Therefore any guaranty given by WPR is limited to the relative reactivity in a standard performed DAS-ELISA with a 10 times diluted positive control originating from Prime Diagnostics.
- By placing any order with WPR, the client declares that he/she has read and accepts these “Special terms of delivery” as well as the “General Terms and Conditions” of WPR and that he/she has complete knowledge of the current product information.
- Client shall use the material in appropriate containment conditions only for research purposes. In no event WPR shall be liable for any use by client of the material or any loss, claim, damage, or liability of whatever kind of nature, which may arise from or in connection with this agreement or the use, handling, storage or transport of the material. WPR shall be safeguarded by the recipient company against any claim regarding these matters. Any damage or loss to the material during transport is at the purchaser's full risk.
- In the event of any liability of WPR, this liability will be limited to the invoice amount for that part of the order to which the liability pertains, on the understanding that this amount shall in no event exceed the amount WPR in such case will receive from its liability insurance.
- Most of the conjugated antibodies are stabilized with bovine serum albumin (BSA). Due to European legislation (EU DIR 1774/2002 and 668/2004) the use of BSA is not allowed in fields of human or veterinary medicine, agriculture, food or cosmetics. By ordering antibodies from WPR client acknowledges that the ordered product, containing BSA, will be used exclusively for research and analytical purposes. If required a form for the written declaration can be obtained on request or downloaded from our website [www.primediagnosics.com](http://www.primediagnosics.com).
- Most products are usually available from stock and are shipped within two weeks on receipt of the order. However, to ensure timely delivery, orders should be placed 8 weeks prior to the desired delivery date. In the event of force majeure WPR will be entitled to suspend performance of the agreement or to terminate the agreement without recourse to the courts and without any liability towards the client. Force majeure on the part of WPR means any circumstance beyond the control of WPR, for example: strikes, fire, war, damage, transport difficulties, export obstructions, defaults on the part of suppliers and legal bars to manufacture or to supply the products.
- Payment should be made by wire transfer (Euro Base Payment applying code SHA (shared costs)) using IBAN of WPR and BIC/SWIFT of RABOBANK within 30 days after date of the invoice.
- Permits, import charges or client's formalities necessary for import by the client of the ordered products are not the responsibility of WPR. It is the responsibility of the client that the client takes care of this before placing an order. All necessary permits need to be in the possession of WPR before an order can be shipped and WPR cannot be held responsible for any damages or import charges that might occur from import problems. The client has to verify if the use of the products is allowed by its government.
- Prices are CPT (Incoterms 2020) in Euros (€), excluding VAT and all orders are charged with € 25.00 for handling. Orders exceeding € 5000.00 are free of handling costs.
- Specificity of the antisera is tested against strains known to be pathogenic at the time of the last actual testing. WPR cannot be held responsible for possible false positive or false negative results caused by newly emerging pathogens, developed resistance, pathogen strains or plant or matrix substances.
- The client has to verify the suitability of the products before purchasing the products. When the product information isn't clearly enough the client has to ask WPR for more information.
- For payments the following VAT number of Wageningen Plant Research have to be used:
  - Domestic: NL811383696B07
  - Foreignn : NL806511618B01.
- WPR does not accept the responsibility for any direct or indirect damage that might arise from the use of delivered products. WPR shall be safeguarded by the recipient company against any claim regarding the delivery and / or the use of the delivered products.


## Schematic overview of the ELISA procedure

These figures do not represent the correct dimensions of the reagents and consumables.



## Schematic overview of the Luminex procedure

These figures do not represent the correct dimensions of the reagents and consumables. The bead mixtures used contain hundreds of beads of different bead addresses. Each bead address is conjugated to one (1) specific antibody for one (1) specific plant pathogenic organisms. This makes it possible to test for multiple pathogens simultaneously, each bead address - plant pathogen combination is specific for the pathogen it was designed for.



Sample	Plant pathogen				
	1 (bead 12)	2 (bead 23)	3 (bead 41)	4 (bead 71)	5 (bead 80)
1	2500	200	1356	230	28
2	27	30	29	33	31
n	x	x	x	x	x


## Schematic overview of the TaqMan procedure

These figures do not represent the correct dimensions of the reagents and consumables.

After DNA isolation specific primers and probe are added to the sample. The probe consists of quencher (Q) and a fluorescent (F) molecule. When fluorophore and quencher are in proximity, quenching inhibits any fluorescence signals.



The above steps can be repeated up to 40 times and in each repetition more fluorescent molecules will be released. In this way the signal is amplified exponentially. When plotted in a graph (Fig.



Comparison between ELISA and Luminex xMAP	
ELISA	Luminex xMAP
<b>Workflow</b>	
1. Coat microtitre plate	1. Sample extraction
2. Sample extraction	2. Add bead mixture to the wells of microtitre plate
3. Transfer samples to the wells of microtitre plate	3. Transfer samples to the wells of microtitre plate
4. Incubate and wash	4. Incubate and wash
5. Add secondary antibody-AP conjugate	5. Add secondary antibody mixture
6. Incubate and wash	6. Incubate and wash
7. Add reporter	7. Add streptavidin-R-phycoerythrin
8. Incubate and analyse	8. Incubate, wash and analyse
<b>Multiplex</b>	
No, multiplex not possible	Yes, up to 20 targets simultaneous detectable
<b>Speed</b>	
Depending on incubation times 16 to 24 hours	Finished in less than 4 hours
<b>Consumables and labour input</b>	
For each new target extra amounts of consumables and labour needed	Consumables and labour input not influenced by number of targets (pathogens)
<b>Flexible</b>	
No, coated or pre-coated microtitre plates needed	Yes, on demand the tests can be performed
Generally, whole plate has to be used	Depending on sample numbers (part of plate can be used)
Evaporation of samples in the wells at the edges of the plate	All wells of microtitre plate can be used
<b>Reliability</b>	
Duplos needed for reliable results	No duplos needed, for each targets more than 50 repetitions are measured


To access the web site of Prime Diagnostics scan the QR code:




Or copy <https://shop.wur.nl/primediagnosics> to your web browser

This will open the webshop of Prime Diagnostics.  
If you already created an account: login

If you are a first time user: Select "Create an Account", and follow the steps in the guide shown below




This will open the following screen:



Fill in the required fields (marked with \*), next click "Create an Account".

After the account is created, please fill in the billing address, if the shipping address is different you need to provide this information as well.



Add address details, and if applicable Chamber of Commerce number

Save address, and go to shipping address

Save shipping address

Once this is completed you will be informed by e-mail that you have created your account.

Before the first order can be made and to confirm that the account was made by you the account needs to be activated. Mail a reply to Prime Diagnostics (primediagnosics@wur.nl) to confirm that your account needs to be activated.

After activation you have full access to the webshop of Prime Diagnostics.

# primediagnos<sup>t</sup>ics

Prime Diagnostics provides high quality reagents for the detection of plant pathogens in various formats:

- ELISA
- Immunofluorescence
- Luminex xMAP and xTAG
- Lateral Flow Devices
- TaqMan and LAMP
- Inocula

From starters to experienced users in the field of plant diagnostics we offer 'state of the art' training with emphasis on improvement of the overall quality of the existing laboratory.

- Companies who want to start up plant diagnostic laboratory activities, Prime Support, a service from Prime Diagnostics can assist in the first steps to start-up plant laboratory activities.
- To existing plant health diagnostics laboratory operations Prime Support offers training to improve skills in plant diagnostics.

Droevendaalsesteeg 1  
P.O. Box 16  
6700 AA Wageningen  
The Netherlands

phone : +31 317 480 613  
e-mail : [primediagnos<sup>t</sup>ics@wur.nl](mailto:primediagnos<sup>t</sup>ics@wur.nl)  
url : [www.primediagnos<sup>t</sup>ics.com](http://www.primediagnos<sup>t</sup>ics.com)  
webshop: [shop.wur.nl/primediagnos<sup>t</sup>ics](http://shop.wur.nl/primediagnos<sup>t</sup>ics)