

FAQ DAS-ELISA

1. *Can I buy a plate reader from PrimeDiagnostics?*

No. You can buy the reagents for detection of plant pathogens from Prime Diagnostics. The machine to measure your microplates can be bought from any supplier of labware.

2. *Do we need a plate washer?*

A plate washer can be used but is not strictly necessary. It is also possible to wash your plates by hand. Hand washes can be done by dipping the plate in demi water or washing buffer or rinsing with tap water or demi water. In case of hand washing take care to empty the plates sufficient.

3. *How do I store the reagents*

Coating, AP-conjugate and VPC's can be stored in the refrigerator. PC's should be stored in a freezer at -20°C. Take care not to use a self-defrosting freezer for storage of any reagent. The temperature changes caused by the defrosting system will harm the reagents.

4. *Do you sell ELISA kits?*

No. We only sell the antibodies and controls for detection of a plant pathogen. Microplates, buffers and substrate are not included.

5. *What kind of antisera do you sell?*

All our antibodies, coatings and AP-conjugates, are polyclonal antibodies raised in rabbit.